Safety, Regulatory and Technical Specification

NeoScan™ 2000

Customer support

For support, please visit our webpage:
neoss.com/neoscan#support

Email support

DE: neoscan.support.de@neoss.com

DK: neoscan.support.dk@neoss.com

ES: neoscan.support.es@neoss.com

FR: neoscan.support.fr@neoss.com

IT: neoscan.support.it@neoss.com

NL: neoscan.support.nl@neoss.com

NO: neoscan.support.no@neoss.com

SE: neoscan.support.se@neoss.com

UK: neoscan.support.uk@neoss.com

US: neoscan.support.us@neoss.com

Other: neoscan.support@neoss.com

NeoScan™ 2000

Notice

The Safety, Regulatory and Technical Specifications User Guide includes information on the safety instructions, regulatory information, and the technical specifications of the devices. We recommend that you thoroughly familiarize yourself with this guide to make the most effective use of your system.

The information in this document is subject to change. Neither Neoss AB nor Neoss Ltd or any of its subsidiaries shall be liable for errors contained herein or for incidental damages in conjunction with the furnishing, performance, or use of this material. No part of this publication may be reproduced without the permission of the manufacturer.

This document is originally written in English.

All trademarks and registered trademarks are the property of their respective holders.

The NeoScan™ 2000 is intended for professional use only. U.S. Federal law restricts this device to sale by or on the order of a dentist.

If any serious incident occurs in relation to the device, the user must report it to Neoss AB and to the competent authority of its Member State in the European Union.

Manual Name: NeoScan 2000 User Manual:
Safety, Regulatory and Technical Specification
Part Number: 16002

Revision: 1

Date: 2023-09

1 – Safety information	5
Indication for use	5
Contraindications	5
Clinical benefits and performance characteristics	5
Warnings and safety instructions	5
Cleaning, disinfecting, sterilizing	8
Clean and disinfect the scanner	8
Clean and sterilize the scanner tip	9
Precautions before use	10
2 – Regulatory information	11
Labels and symbols	11
Regulatory information	12
Compliance with European and international standards	12
Classification in accordance with EN/IEC 60601-1	13
Conformity with EN/IEC 60601-1-2	13
Guidance and manufacturer's declarations	15
Accessories	16
Other equipment	16
3 – Technical specification	17
Model	17
Technical specifications	17
Environmental requirements	17

1 – Safety information

Indication for use

The NeoScan 2000 Intraoral Scanner (herein after referred as 'scanner') is a digital optical scanning device used to obtain digital impressions of hard and soft tissues such as teeth, gums, and mucous membrane by means of oral scanning, for oral restoration and orthodontic treatment of malocclusion.

The scanner could be used for both adults and children in clinical practice.

Contraindications

None

Clinical benefits and performance characteristics

The scanner benefits dental practices by enabling practitioners to acquire digital impressions with the quality and accuracy required for digital CAD/CAM dental applications. The actual performance of the device is dependent on the user's training and operating execution. The user is solely responsible for the accuracy, completeness, and adequacy of the acquired data.

The scanner is designed to acquire 3D models in the followings:

- Upper jaw
- Lower jaw
- Buccal bite registration

Warnings and safety instructions

DANGER OF ELECTRIC SHOCK

This is an electrical unit. Do **NOT** expose it to water spray. Such action can cause an electric shock or a malfunction of the unit.

CAUTION: All known residual risks, contraindications, or undesirable side effects are listed in this document. If any serious incident occurs in relation to the device, you must report it to Neoss AB and to the competent authority of your Member State in the European Union.

WARNING

Scanner

- You **MUST** read and understand this safety information before using the scanner.
- This scanner shall only be used inside hospitals and other professional healthcare facilities and **MUST NOT** be used near high frequency surgical equipment and the RF shielded room of an ME System for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.
- Before using the scanner, check the outer surfaces of the unit and any accessories to ensure there are no rough surfaces, sharp edges, or protrusions which may cause a safety hazard.
- You are responsible for the operation and maintenance of the scanner. You **MUST** have training to use the scanner.
- **DO NOT** place objects within the field of operation of the unit.
- When the unit is not in use, ensure that the scanner is turned OFF.
- **DO NOT** use the scanner in conjunction with oxygen-rich environments. This unit is not intended for use with flammable anesthetics or flammable agents.

- **DO NOT** pull or twist the cable.
- **DO NOT** drop the scanner or the accessories.
- **DO NOT** sterilize the scanner.
- **DO NOT** expose the scanner to water spray or submerge it in water or disinfectant.
- **DO NOT** expose the scanner to high vibrations.
- **DO NOT** expose the scanner to ultraviolet radiation directly. The scanner is not designed for ultraviolet disinfection.
- **DO NOT** stare at the LED emission window.
- When the Scanner Tip is removed, install the front protective cover to protect the Lens Window.
- **DO NOT** remove the cover of any scanner components. The scanner contains no user-serviceable parts. For any repairs, contact a qualified Neoss AB service technician.
- **DO NOT** replace the cables provided with the scanner with other cables. Doing so may damage the scanner and adversely affect the safety protection and EMC performance of the scanner.
- Any other equipment not complying with IEC-60601 shall be kept at least 1.5 meters away from the patient.
- If the equipment is faulty, turn it OFF, display an "Out of Service" notice, and contact a qualified Neoss AB service technician.
- Using components, accessories, cables, and spare parts other than those specified or provided by the manufacturer of this equipment may impair the safety protection of the scanner and may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- No modification of this equipment is allowed.
- Additional multiple outlet strips or extension cords should not be connected to the system.
- The maximum temperature of the applied part may reach to 43°C; to avoid overheating, do not use it for extended periods.
- To isolate the holder from power supply, unplug the USB connector from the USB port.
- **DO NOT** maintain or service this equipment while it is in use with the patient.
- Connection of the PEMS (Programmable Electrical Medical System) to an IT NETWORK that includes other equipment could result in risks to patients, operators, or third parties. The responsible organization should identify, analyze, evaluate, and control these risks.
- Patients with oral mucosal disease, mental illness, severe respiratory disease, asthma, Parkinson's disease, hyperactivity disease are forbidden.
- Patients with moderate or severe opening limitation should use it with caution.

Computer

- **DO NOT** place any equipment which does not comply with IEC 60601-1 in the immediate vicinity of the patient. Leave at least 1.5 meters distance between the patient and the equipment.
- The scanner is only intended to be connected to a computer that is at least IEC 60950/IEC 62368, or equivalent standards certified. Connecting the scanner to other equipment may be hazardous.
- See the installation guide for your computer for information about the data processing system, computer, and screen. Leave a sufficient amount of clear space around the computer to ensure that it is properly ventilated.
- Position the screen to avoid light reflections from internal or external lighting for maximum image quality and visual comfort.

Battery

- Do not dismantle, open, or shred secondary cells or batteries.
- Do not expose cells or batteries to heat or fire. Avoid storage in direct sunlight.
- Do not short-circuit a cell or a battery. Do not store cells or batteries haphazardly in a box or drawer where they may short-circuit each other or be short-circuited by other metal objects.
- Do not remove a cell or battery from its original packaging until required for use.
- Do not subject cells or batteries to mechanical shock.
- In the event of a cell leaking, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice.
- Do not use any charger other than that specifically provided for use with the equipment. Refer to the manufacturer's instructions for proper charging instructions.
- Do not use any cell or battery which is not designed for use with the equipment.
- Always purchase the battery recommended by the device manufacturer for the equipment.
- Keep cells and batteries clean and dry.
- Wipe the cell or battery terminals with a clean, dry cloth if they become dirty.
- Do not leave a battery on prolonged charge when not in use.
- After extended periods of storage, it may be necessary to charge and discharge the cells or batteries several times to obtain maximum performance.
- Retain the original product literature for future reference.
- Use only the cell or battery in the application for which it was intended.
- When possible, remove the battery from the equipment when not in use.
- Please remove battery using specified tools. Removing the battery by hand is not allowed.
- Dispose of the battery properly.

Disposal

This equipment contains certain materials and chemical compounds incidental to the manufacture of electrical and electronic equipment, and improper "end-of-life" disposal of such equipment can result in environmental contamination. Therefore, this equipment should not be disposed of as ordinary household waste but should instead be delivered to a designated electrical and electronic waste disposal or recycling center. For further information on disposing of electrical and electronic waste, contact the cognizant authority within the local jurisdiction.

Dispose of the Scanner Tips according to standard operating procedures or local regulations for the disposal of contaminated medical waste. For additional Scanner Tips, contact your dealer.

Cleaning, disinfecting, sterilizing

Perform the following maintenance activities on your scanner and accessories regularly.

To ensure maximum hygienic safety for the patient, carefully follow the instructions to prepare the scanner for use.

To ensure maximum hygienic safety for the patient and to minimize the risk of cross-contamination, carefully perform the following maintenance activities on your scanner and accessories.

After each patient:

- Clean and disinfect the scanner. See “Clean and Disinfect the Scanner”.
- Clean the Scanner Tip, and then perform Autoclave Sterilization (See “Clean the Scanner Tip” and “Sterilize the Scanner Tip”).

There are 2 models of Scanner Tip:

Model	Size	UDI-DI	Manual Cleaning	134° C Sterilization
TP102N	Large	(01)06973993441034	✓	✓
TP103N	Small	(01)06973993441041	✓	✓

Clean and disinfect the scanner

General Warnings

- Read and follow the warnings and personal protection instructions provided in the Safety Data Sheet (SDS) for the disinfectant used to process the scanner.
- You must wear gloves while cleaning and disinfecting the scanner.
- The scanner must be disinfected with a recommended intermediate-level disinfectant solution with tuberculocidal activity between patients.
- **DO NOT** use a disinfectant containing phenolics or iodophors; doing so will damage the surface coating of the scanner.
- Never put the scanner in an autoclave device or immerse it in water or the disinfectant solution.
- Excessive fluids can damage the scanner.
- Do not use cotton, cloth, or tissues soaked with disinfectant to disinfect the scanner.

Clean the Scanner

- If the scanner is visibly contaminated with blood and/or body fluids, you must clean it before disinfecting it.
- To clean the scanner, follow these steps:
 - Dampen (do not soak) a lint-free cloth with lukewarm water.
 - Remove the blood and/or body fluids with the dampened lint-free cloth.

Disinfect the Scanner

After each patient, the scanner must be thoroughly disinfected.

To adequately disinfect the scanner, follow the disinfectant manufacturer's instructions for the appropriate contact time.

CAUTION: If the scanner is visibly soiled, it must be thoroughly cleaned prior to disinfecting. See "Clean the Scanner"

To disinfect the scanner, follow these steps:

- Remove the Scanner Tip.
- Remove all visible soil (see "Clean the Scanner").
- Use a commercially prepared intermediate level disinfectant wipe. Follow the manufacturer's instructions for contact time.

Recommended disinfectant wipes: CaviWipes™*

WARNING: Using a disinfectant that has not been approved may cause damage to the scanner.

Thoroughly wipe all surfaces of the scanner, **DO NOT** allow liquid to enter through the gap, air outlet, or pin holes.

WARNING: Do not rinse.

- Allow to air dry.
- After the scanner has dried, use a clean, lint-free cloth dampened with water to remove residual disinfectant from the surface of the scanner.

Clean and sterilize the scanner tip

WARNING

- Wear gloves when handling a contaminated Scanner Tip.
- Read and follow the warnings and personal protection instructions provided in the manufacturer's SDS for the detergent used to clean the Scanner Tip.
- Do not soak the Scanner Tip in disinfectant for a long period.
- Dry the Scanner Tip thoroughly before mounting onto the scanner.
- Do not use an ultrasonic cleaning machine to clean the Scanner Tip.
- Do not soak the Scanner Tip with alcohol-based disinfectants.

Clean the Scanner Tip

To manually cleaning the Scanner Tip, follow these steps:

- Rinse excess soil from the Scanner Tip (2 minutes).
- Using a soft brush, apply enzymatic detergent solution (e.g., Metrex EmPower) to all surfaces.
- Rinse under clean, running water (2 minutes).
- Inspect Scanner Tip. If the Scanner Tip is not clean, repeat the last three steps.
- Use a lens tissue or lint-free cloth to remove any dust from the mirror in the Scanner Tip.

* CaviWipes are a trademark of Metrex Research, LLC.

Sterilize the Scanner Tip

Scanner Tips shipped non-sterile. You must sterilize them before use.

CAUTION: If you limit the exposure time at 134°C to not more than 6 minutes, you can autoclave the Scanner Tip up to 60 cycles.

To sterilize the Scanner Tip, follow these steps:

- Put the Scanner Tip into a sealed steam sterilization pouch.
- Place the Scanner Tip in a steam autoclave for sterilization.
 - For TP102N and TP103N, exposure temperature should be set as 134°C.
 - For TP102N and TP103N, exposure time should exceed 3 minutes.
 - For TP102N and TP103N, exposure time should not exceed 6 minutes.

Precautions before use

Perform the following activities on your product and accessories before use.

Visually Inspecting the Scanner for Damage

Visually inspect the scanner for damage or signs of deterioration by doing the following:

- Inspect the Lens Window.
- Inspect around the buttons and contact points.

If damage is noted, do not use the scanner and contact your representative.

Visually Inspecting the Scanner Tip

Visually inspect the Scanner Tip for signs of deterioration by doing the following:

- Verify that the Scanner Tip is not damaged, and its components are not detached.
- Verify that the Scanner Tip mirror does not have any smudges or scratches on it.

If deterioration is noted, replace the Scanner Tip.

WARNING

- The Lens Window on the scanner is a delicate optical component. Mount the front Protective Cover to protect the Lens Window from damage and dirt when the scanner is not in use.
- The mirror in the Scanner Tip is a delicate optical component. Its clean and undamaged surface is critical to scan quality.

In the event that you see poor scan quality or an unclear video preview in the software, clean the Scanner Tip mirror and the Lens Window of the Handpiece using a microfiber cleaning swab, applying ethanol that is free of impurities.

2 – Regulatory information

Labels and symbols

	Type BF applied part
	Separate collection for waste electric and electronic equipment (WEEE) is required.
	Manufacturer
	Date of manufacture
	Importer
	Refer to instruction manual/booklet
	Direct current
	Sterilizable in a steam sterilizer (autoclave) at 134°C
	Model number
	Serial number
	Rechargeable battery
	CE marking applicable for European Union
	This device complies with part 15 of the FCC rules
	Authorized Representative in the European Community

Regulatory information

The device complies with the following regulations:

- MDR: (EU) 2017/745 Medical Device Regulation, Class I following the Rule 5.
- FDA Center for Devices & Radiological Health CDRH – Title 21 CFR 872.3661 (USA).
- Medical Devices Regulations (Canada).
- RoHS: Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment, 2011/65/EU Annex II and its amendment Directive (EU) 2015/863.
- RED: Directive 2014/53/EU The Radio Equipment Directive.
- FCC: Part 15 of The Federal Communications Commission Rules.
- ISED: Innovation, Science and Economic Development Canada.

Compliance with European and international standards

EN/IEC 60601-1:	Medical Electrical Equipment, Part 1: General Requirements for Basic Safety and Essential Performance
ANSI/AAMI ES 60601-1:	Medical Electrical Equipment – Part 1: General requirements for basic safety and essential performance
CAN/CSA-C22.2 NO. 60601-1:	Medical Electrical Equipment – Part 1: General requirements for basic safety and essential performance
EN/IEC 60601-1-2:	Medical Electrical Equipment, Part 1–2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests
EN/IEC 80601-2-60:	Medical electrical equipment – Part 2–60: Particular requirements for the basic safety and essential performance of dental equipment
EN/IEC 62471:	Photobiological safety of lamps and lamp systems
EN/ISO 17664:	Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices
EN/ISO 17665-1:	Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices
EN/IEC 60601-1-6:	Medical Electrical Equipment, Part 1–6: General requirements for basic safety and essential performance – Collateral Standard: Usability
EN/IEC 62366-1:	Medical devices – Part 1: Application of usability engineering to medical devices
EN/IEC 62304:	Medical device software – Software life cycle Processes
EN ISO 10993:	Biological evaluation of medical devices
ISO 14971:	Medical devices – Application of risk management to medical devices
EN/ISO 15223-1:	Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements

EN ISO 20417:	Medical devices – Information to be supplied by the manufacturer
ISO 9687:	Dentistry – Graphical symbols for dental equipment
AAMI TIR 12:	Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers
AAMI TIR 30:	A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices
EN/IEC 62133-2:	Secondary cells and batteries containing alkaline or other non-acid electrolytes – Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications – Part 2: Lithium systems
EN 50566:	Product standard to demonstrate the compliance of wireless communication devices with the basic restrictions and exposure limit values related to human exposure to electromagnetic fields in the frequency range from 30 MHz to 6 GHz: hand-held and body mounted devices in close proximity to the human body.
EN 301489-1:	Electromagnetic compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements; Harmonized Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU and the essential requirements of article 6 of Directive 2014/30/EU
EN 301489-17:	Electromagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems; Harmonized Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU
EN 301893:	5 GHz RLAN; Harmonized Standard covering the essential requirements of article 3.2 of Directive 2014/53/EU

Classification in accordance with EN/IEC 60601-1

Type of protection against electric shock: Internally powered

Degree of protection against electric shock: Type BF Applied Part

Mode of operation: Continuous operation

Flammable anesthetics: Not suitable for use in the presence of flammable anesthetics or a mixture of flammable anesthetics with air or oxygen or nitrous oxide.

Conformity with EN/IEC 60601-1-2

IEC 60601-1-2 EMC requirements and tests, Medical Electrical Equipment including CISPR 11 Group 1, Class B.

Electromagnetic Compatibility Precautions

Medical electrical equipment requires special precautions regarding electromagnetic compatibility (EMC). Medical equipment must be installed and put into service according to the EMC information provided in this documentation.

Other equipment can interfere with communications with the device, even if the equipment complies with CISPR emissions requirements.

WARNING: Portable and mobile RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the DEVICE, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could happen.

WiFi

The device operates with the 802.11a/n/ac protocol. Channel 38 or 46 is priority used for the handpiece. The channel bandwidth is 40 MHz. The frequency range is 5150–5250 MHz, 5725–5850 MHz (the actual frequencies are dependent on local regulations and the configuration of the product). The maximum output power is 17.88 dBm.

This device complies with part 15 of the FCC Rules and contains license exempt transmitter(s)/ receiver(s) that comply with ISSED's license-exempt RSS(s).

Operation is subject to the following two conditions:

- This device may not cause interference.
- This device must accept any interference, including interference that may cause undesired operation of the device.

NOTE: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

SAR Value for Handpiece:

0.798 W/Kg, 10 g for CE

0.011 W/Kg, 1 g for FCC

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Limited by local law regulations, the version for North America does not have a region selection option.

The device is for indoor use only and operates in the 5150–5250 MHz band to reduce the potential for harmful interference to co-channel mobile satellite systems.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not in-stalled and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Restrictions in the 5 GHz band: According to Article 10 (10) of Directive 2014/53/EU, the packaging shows that this radio equipment will be subject to some restrictions when placed on the market in Belgium (BE), Bulgaria (BG), the Czech Republic (CZ), Denmark (DK), Germany (DE), Estonia (EE), Ireland (IE), Greece (EL), Spain (ES), France (FR), Croatia (HR), Italy (IT), Cyprus (CY), Latvia (LV), Lithuania (LT), Luxembourg (LU), Hungary (HU), Malta (MT), Netherlands (NL), Austria (AT), Poland (PL), Portugal (PT), Romania (RO), Slovenia (SI), Slovakia (SK), Finland (FI), Sweden (SE), Turkey (TR), Norway (NO), Switzerland (CH), Iceland (IS), and Liechtenstein (LI). The WLAN function for this device is restricted to indoor use only when operating in the 5150 to 5250 MHz frequency range.

AT, BE, BG, CH, CY, CZ, DE, DK, EE, EL, ES, FI, FR, HR, HU, IE, IS, IT, LI, LT, LU, LV, MT, NL, NO, PL, PT, RO, SE, SI, SK, TR

Guidance and manufacturer's declaration

Guidance and Manufacturer's Declaration – Electromagnetic Emission

The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1 Class B	The DEVICE uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity for Equipment and Systems (from Non-RF wireless communication equipment)

The DEVICE is intended for use in the electromagnetic environment specified below (IEC 60601-1-2 Table 4). The customer or the user of the DEVICE should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz	3 V/m 80 MHz – 2.7 GHz	Environment of a professional healthcare facility.

NOTE: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DEVICE is used exceeds the applicable RF compliance level above, the DEVICE should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the DEVICE.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity for Equipment and Systems (from RF wireless communication equipment)

For the immunity to proximity fields from RF wireless communications equipment, the DEVICE is compliant with the test levels specified below, according to IEC60601-1-2 Table 9. The customer or user of the DEVICE should assure that it is used in such an environment.

Test Frequency (MHz)	Band (MHz)	Immunity Test Levels
385	380 – 390	Pulse modulation 18 Hz, 27 V/m
450	430 – 470	FM, ±5 kHz deviation, 1 kHz sine, 28 V/m
710	704 – 787	Pulse modulation 217 Hz, 9 V/m
745		
780		
810	800 – 960	Pulse modulation 18 Hz, 28 V/m
870		
930		
1720	1700 – 1990	Pulse modulation 217 Hz, 28 V/m
1845		
1970		
2450	2400 – 2570	Pulse modulation 217 Hz, 28 V/m
5240	5100 – 5800	Pulse modulation 217 Hz, 9 V/m
5500		
5785		

Accessories

The use of accessories other than those specified, except for those sold by the manufacturer of the equipment, as replacement parts for internal components may result in increased emissions or decreased immunity of the medical equipment.

Other equipment

 **WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify normal operation.

3 – Technical specification

Model

NeoScan 2000

Technical specifications

Components	Technical Specifications
Weight	Scanner (including battery, without Scanner Tip): 245 g Tip Large – Sterilization (TP102N): 13 g Tip Small – Sterilization (TP103N): 12 g
Color	3D full color
USB wireless Card Port	USB 3.0
USB charge Port	Type C Port
Field of view	Tip Large – Sterilization (TP102N): 16 mm × 14 mm Tip Small – Sterilization (TP103N): 12 mm × 12 mm
Depth of view	15mm
Configuration requirement of workstation	Processor: Intel® Core™ i7 9th Generation, base frequency 2.6 GHz (or better) Memory: 32 GB (or more) DDR4, frequency 2666 MHz (or better) Disk: 512G (or more) SSD Graphics card: NVIDIA® GeForce® GTX 1650 (or better) Display: 15.6" FHD (1920 × 1080) (or more) Others: USB 3.0 port Operating system: Windows 10 (build 18362 or newer)/ Win 11, 64 bit Optional: Touch screen

 **CAUTION:** It is MANDATORY to check that your system configuration is compatible with the computer system requirements for the accompanying software.

Environmental requirements

Components	Environmental Requirements
Operating temperature	15°C ~ 30°C
Transport and storage temperature	-10°C ~ 50°C
Operating relative humidity	10% ~ 65% RH
Transportation and storage relative humidity	10% ~ 95% RH
Operating atmospheric pressure	70 ~ 106 KPa
Transportation and storage atmospheric pressure	60 ~ 106 KPa

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Distributed by:

Neoss Ltd. and its subsidiaries as above

Manufacturer:

Alliedstar Medical Equipment Co., Ltd. No.1-4, Floor 3, Unit 2, Building 3, No. 222, West third section, Waihuan Road, Yanjiang District, Ziyang, 641300, Sichuan, P.R. China | Email: support@allied-star.com

EU Authorized Representatives:

SUNGO Europe BV., Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands

Neoss AG

The Circle 23, 8058 Zurich-Airport, Switzerland

E: distributors@neoss.com

E: info@neoss.com

For support, please visit our webpage:

neoss.com/neoscan

NeoScan customer support: +44 1423 533 991

Contact your local representative for product availability details. All products are not available in every market.

16002_1 EN, DA, FR, DE, IT, NO, PT, ES, SV 2023-09 © Neoss Group, 2023. Copyrights, design rights and trademarks, Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorization of Neoss Group. Neoss and the N logo are registered trademarks of Neoss Group. NeoScan™, NeoPro™ and NeoConnect™ are trademarks of Neoss Group. CaviWipes are a trademark of Metrex Research, LLC.

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Neoss AB, Arvid Wallgrens backe 20, SE-413 46 Göteborg | Phone: +46 (0)31 88 12 80 | Email: info@neoss.se