



Product Catalogue

2022 | 2023




Intelligent Simplicity

Explore Intelligent Simplicity

Founded in 2000, Neoss® offers intelligent solutions that are intuitively easy to use. Our products allow dental professionals to provide reliable and cost-effective treatments to their patients with predictable long-term results. We strive to set new standards by leading the market with ingenuity and integrity. In developing smart treatment solutions and working closely with each practice, Neoss makes the complex less complicated. We call that "Intelligent Simplicity".

Our NeossONE® solution; one prosthetic platform across three implant ranges, means that you can use the same abutment for 70 different implants, providing you with a more cost-effective solution, and it makes training the whole dental team easier. This is the essence of our Intelligent Simplicity ethos.

We believe it's part of our remit to help dental teams explain to patients why dental implants are the optimum aesthetic and practical solution and what makes the Neoss system a great investment. In developing smarter treatments and working closely with you, Neoss is the ultimate practice-builder partner.



Implants

- 4 Neoss ProActive® Edge Implants
- 5 Neoss ProActive® Straight Implants
- 6 Neoss ProActive® Tapered Implants
- Neoss ProActive® Wide Implants
- Neoss ProActive® Sinus Implants
- 7 Neoss ProActive® Straight NP Ø3.25 Implants

Biomaterials

- 8 **Bone Graft Materials**
 - THE Graft™ in vials
 - THE Graft™ in syringes
 - THE Graft™ in blocks
- Resorbable Membranes**
 - Collagen Membranes
 - Wound Dressings
- 9 **Non-Resorbable Membranes**
 - Ti-reinforced PTFE Membranes
 - Non-reinforced PTFE Membranes
- Fixation Kit**
- Sutures**
 - Biotex®
 - Other

Digital

- 10 **NeoScan**
 - Digital Impression**
 - Esthetic Healing Abutments with ScanPeg
 - ScanPeg
 - Scan Body
- 11 **Analogs**
 - Digital Prosthetics**
 - Provisional Abutments
 - NeoBase® Abutments SSC (Straight Screw Channel)
- 12 NeoBase® Abutments ASC (Angulated Screw Channel)
- Neoss TiBase Abutments (compatible with Dentsply Sirona®)
- Custom Solutions**
 - Neoss NT Blanks
- 13 Neoss Custom Prosthetics – powered by ARC

Prosthetics

- 14 **Healing Abutments**
 - Esthetic Healing Abutments with ScanPeg
 - Healing Abutments
- Provisional Abutments**
 - Provisional Titanium Abutments
- 15 Esthetic Tissue Formers
- Final Abutments**
 - Prepable Titanium Abutments
 - Prepable Titanium Abutments, Alternative Emergence Profiles
- 16 Zirconia Abutments
- Final Restorations – Casting**
 - Gold NeoLink® (Direct Casting)
 - CoCr Abutments (Indirect Casting)
 - Burnout Abutments (Direct Casting)
- 17 Ti NeoLink®
- Multiple Unit Abutments**
 - Access Abutments (Screw-retained multiple unit restoration)
- 18 **Overdenture Abutments**
 - Equator Abutments
- 19 Ball Abutments
- Bar Abutment
- Impressions & Replicas**
 - Impressions Copings
 - Replicas
- 20 **Screws**

Instruments & More

- 21 **Drills**
- 22 **Kits**
 - Surgical and Prosthetic Trays
- 23 Trays & Instrument Kits
 - Guide
 - Neoss® Guide Kit – SPARE PARTS
 - Demo Kit
- 24 **Instruments**
 - Surgical Equipment**
 - Neoss Penguin^{®FA}

Support

- 25 Prosthetic Solution Standard and Narrow Platform, overview
- Torque and Speed recommendation, overview
- 26 Symbols, names and abbreviations
- You can now buy your products online!
- Disclaimer of Liability
- 28 Contact details

Neoss ProActive® surface – *Proven predictability* in implant treatment for all patients since 2009



ProActive is a modern super-hydrophilic and ultra-clean implant surface, that was introduced already in 2009. ProActive has a proven clinical success and is backed-up by long-term data¹.

✓ Ultra-clean, super-hydrophilic and magnesium-ion enriched


- Accelerated and increased strength of osseointegration²
- Maintained implant stability during healing³

✓ Defined dual surface (blasted and etched) with lower roughness at implant collar

- Stable bone levels over time – Bone maintenance⁴
- 10-year clinical data on > 1200 patients showing > 98% cumulative survival rate^{4,5,6}




1) Engman F et al. *Letters Impl Dent.* 2022;2:6-11. 2) Sennerby L et al. *Letters Impl Dent.* 2017;1:7-10.
3) Vanden Bogaerde L et al. *Int J Dent.* 2016;84:24931. 4) Sahlin H et al. *Letters Impl Dent.* 2022;2:18-24.
5) Sennerby L et al. *Letters Impl Dent.* 2022;2:29-32. 6) Hassfurther N et al. *Clin Oral Implants Res.* 2021;32(S22):156-157.



Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments in PEEK in two different sizes and Healing Abutment Screw.





Neoss ProActive® Edge Implants



Art. No.	Description	A	B
21301	Implant Kit – Neoss ProActive® Edge Ø3.5×9mm		
21302	Implant Kit – Neoss ProActive® Edge Ø3.5×11mm	4.0mm	1.6mm
21303	Implant Kit – Neoss ProActive® Edge Ø3.5×13mm		
21304	Implant Kit – Neoss ProActive® Edge Ø3.5×15mm		
21305	Implant Kit – Neoss ProActive® Edge Ø4.0×9mm		
21306	Implant Kit – Neoss ProActive® Edge Ø4.0×11mm	4.3mm	1.9mm
21307	Implant Kit – Neoss ProActive® Edge Ø4.0×13mm		
21308	Implant Kit – Neoss ProActive® Edge Ø4.0×15mm		
21309	Implant Kit – Neoss ProActive® Edge Ø4.5×9mm		
21310	Implant Kit – Neoss ProActive® Edge Ø4.5×11mm	4.9mm	2.2mm
21311	Implant Kit – Neoss ProActive® Edge Ø4.5×13mm		
21312	Implant Kit – Neoss ProActive® Edge Ø4.5×15mm		
21313	Implant Kit – Neoss ProActive® Edge Ø5.0×9mm		
21314	Implant Kit – Neoss ProActive® Edge Ø5.0×11mm	5.4mm	2.4mm
21315	Implant Kit – Neoss ProActive® Edge Ø5.0×13mm		
21316	Implant Kit – Neoss ProActive® Edge Ø5.0×15mm		



SP = Standard Platform
Ø3.5–6.5 Implants



Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments (Ø5.0–6.5 implant kits include a wider version) in PEEK in two different sizes and Healing Abutment Screw.

Neoss ProActive® Straight Implants

Art. No.	Description	A	B
21181	Implant Kit – Neoss ProActive® Straight Ø3.5×7 mm		
21182	Implant Kit – Neoss ProActive® Straight Ø3.5×9 mm		
21183	Implant Kit – Neoss ProActive® Straight Ø3.5×11 mm	4.0 mm	2.4 mm
21184	Implant Kit – Neoss ProActive® Straight Ø3.5×13 mm		
21185	Implant Kit – Neoss ProActive® Straight Ø3.5×15 mm		
21186	Implant Kit – Neoss ProActive® Straight Ø3.5×17 mm		
21187	Implant Kit – Neoss ProActive® Straight Ø4.0×7 mm		
21188	Implant Kit – Neoss ProActive® Straight Ø4.0×9 mm	4.0 mm	2.8 mm
21189	Implant Kit – Neoss ProActive® Straight Ø4.0×11 mm		
21190	Implant Kit – Neoss ProActive® Straight Ø4.0×13 mm		
21191	Implant Kit – Neoss ProActive® Straight Ø4.0×15 mm		
21192	Implant Kit – Neoss ProActive® Straight Ø4.0×17 mm		
21193	Implant Kit – Neoss ProActive® Straight Ø4.5×7 mm		
21194	Implant Kit – Neoss ProActive® Straight Ø4.5×9 mm	4.5 mm	3.2 mm
21195	Implant Kit – Neoss ProActive® Straight Ø4.5×11 mm		
21196	Implant Kit – Neoss ProActive® Straight Ø4.5×13 mm		
21197	Implant Kit – Neoss ProActive® Straight Ø4.5×15 mm		
21198	Implant Kit – Neoss ProActive® Straight Ø4.5×17 mm		
21199	Implant Kit – Neoss ProActive® Straight Ø5.0×7 mm		
21200	Implant Kit – Neoss ProActive® Straight Ø5.0×9 mm	5.0 mm	3.7 mm
21201	Implant Kit – Neoss ProActive® Straight Ø5.0×11 mm		
21202	Implant Kit – Neoss ProActive® Straight Ø5.0×13 mm		
21203	Implant Kit – Neoss ProActive® Straight Ø5.0×15 mm		
21205	Implant Kit – Neoss ProActive® Straight Ø5.5×7 mm		
21206	Implant Kit – Neoss ProActive® Straight Ø5.5×9 mm	5.5 mm	4.2 mm
21207	Implant Kit – Neoss ProActive® Straight Ø5.5×11 mm		
21208	Implant Kit – Neoss ProActive® Straight Ø5.5×13 mm		
21252	Implant Kit – Neoss ProActive® Ø6.0×7 mm	6.0 mm	4.5 mm
21211	Implant Kit – Neoss ProActive® Ø6.5×7 mm	6.7 mm	



ProActive Straight Ø3.5



ProActive Straight Ø4.0



ProActive Straight Ø4.5




ProActive Straight Ø5.0



ProActive Straight Ø5.5



Short Implant



Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments (Ø5.0 – 5.5 implant kits include a wider version) in PEEK in two different sizes and Healing Abutment Screw.



ProActive Tapered Ø3.5



ProActive Tapered Ø4.0



ProActive Tapered Ø4.5



ProActive Tapered Ø5.0



ProActive Tapered Ø5.5



ProActive Ø6.0



ProActive Ø6.5



ProActive Ø6.5



Neoss ProActive® Tapered Implants


Art. No.	Description	A	B
21221	Implant Kit – Neoss ProActive® Tapered Ø3.5×9 mm		
21222	Implant Kit – Neoss ProActive® Tapered Ø3.5×11 mm	4.0 mm	2.1 mm
21223	Implant Kit – Neoss ProActive® Tapered Ø3.5×13 mm		
21224	Implant Kit – Neoss ProActive® Tapered Ø3.5×15 mm		
21227	Implant Kit – Neoss ProActive® Tapered Ø4.0×9 mm		
21228	Implant Kit – Neoss ProActive® Tapered Ø4.0×11 mm	4.3 mm	2.3 mm
21229	Implant Kit – Neoss ProActive® Tapered Ø4.0×13 mm		
21230	Implant Kit – Neoss ProActive® Tapered Ø4.0×15 mm		
21233	Implant Kit – Neoss ProActive® Tapered Ø4.5×9 mm		
21234	Implant Kit – Neoss ProActive® Tapered Ø4.5×11 mm	4.9 mm	2.8 mm
21235	Implant Kit – Neoss ProActive® Tapered Ø4.5×13 mm		
21236	Implant Kit – Neoss ProActive® Tapered Ø4.5×15 mm		
21239	Implant Kit – Neoss ProActive® Tapered Ø5.0×9 mm		
21240	Implant Kit – Neoss ProActive® Tapered Ø5.0×11 mm	5.4 mm	3.3 mm
21241	Implant Kit – Neoss ProActive® Tapered Ø5.0×13 mm		
21242	Implant Kit – Neoss ProActive® Tapered Ø5.0×15 mm		
21245	Implant Kit – Neoss ProActive® Tapered Ø5.5×9 mm		
21246	Implant Kit – Neoss ProActive® Tapered Ø5.5×11 mm	5.9 mm	3.8 mm
21247	Implant Kit – Neoss ProActive® Tapered Ø5.5×13 mm		

Neoss ProActive® Wide Implants

Art. No.	Description	A	B
21252	Implant Kit – Neoss ProActive® Ø6.0×7 mm		4.5 mm
21250	Implant Kit – Neoss ProActive® Ø6.0×9 mm	6.0 mm	4.3 mm
21251	Implant Kit – Neoss ProActive® Ø6.0×11 mm		4.3 mm


Neoss ProActive® Sinus Implants

Art. No.	Description	A
21211	Implant Kit – Neoss ProActive® Ø6.5×7 mm	
21212	Implant Kit – Neoss ProActive® Ø6.5×9 mm	6.7 mm
21213	Implant Kit – Neoss ProActive® Ø6.5×11 mm	



Neoss ProActive® NP Ø3.25 Implant

The Neoss narrow platform implant has a Ø3.25mm thread diameter with a Ø3.5mm restorative platform. Key NP components and labels are colour coded in royal blue.



Neoss Implant Kit

The Neoss ProActive® NP Ø3.25 implant kit includes Implant, Cover Screw, Healing Abutment and Healing Abutment Screw.

Neoss ProActive® Straight NP Ø3.25 Implants

Art. No.	Description	A	B
21176	Implant Kit – Neoss ProActive® Straight, NP Ø3.25 × 9 mm		
21177	Implant Kit – Neoss ProActive® Straight, NP Ø3.25 × 11 mm	3.5mm	2.2mm
21178	Implant Kit – Neoss ProActive® Straight, NP Ø3.25 × 13 mm		
21179	Implant Kit – Neoss ProActive® Straight, NP Ø3.25 × 15 mm		



ProActive Straight NP Ø3.25

Solutions for a variety of *maxillofacial procedures*

Bone Graft Materials

THE Graft™ in vials



THE Graft™ in vials

Art. No.	Description	Weight	Volume
BG-A15	THE Graft Granules 0.25 – 1.00 mm	0.15 g	0.36 cc
BG-A25	THE Graft Granules 0.25 – 1.00 mm	0.25 g	0.60 cc
BG-A05	THE Graft Granules 0.25 – 1.00 mm	0.50 g	1.20 cc
BG-A10	THE Graft Granules 0.25 – 1.00 mm	1.00 g	2.40 cc
BG-A20	THE Graft Granules 0.25 – 1.00 mm	2.00 g	4.80 cc
BG-B05	THE Graft Granules 1.00 – 2.00 mm	0.50 g	1.80 cc
BG-B10	THE Graft Granules 1.00 – 2.00 mm	1.00 g	3.60 cc

THE Graft™ in syringes



THE Graft™ in syringes

Art. No.	Description	Weight	Volume
TG-AS25	THE Graft Syringe 0.25 – 1.00 mm	0.10 g	0.25 cc
TG-AS05	THE Graft Syringe 0.25 – 1.00 mm	0.21 g	0.50 cc
TG-AS10	THE Graft Syringe 0.25 – 1.00 mm	0.42 g	1.00 cc
TG-BS25	THE Graft Syringe 1.00 – 2.00 mm	0.07 g	0.25 cc
TG-BS05	THE Graft Syringe 1.00 – 2.00 mm	0.14 g	0.50 cc
TG-BS10	THE Graft Syringe 1.00 – 2.00 mm	0.28 g	1.00 cc

THE Graft™ in blocks



THE Graft Collagen Block

Art. No.	Description	Weight	Volume
TCB-01	THE Graft Collagen Block	0.15 g	7×7×7 mm
TCB-02	THE Graft Collagen Block	0.3 g	8×9×10 mm
TCB-03	THE Graft Collagen Block	0.6 g	10×11×12 mm

Resorbable Membranes

Collagen Membranes



Collagen Flex

Collagen Firm


Art. No.	Description	A	B
64050	NeoGen® Collagen Membrane Flex S	20 mm	15 mm
64051	NeoGen® Collagen Membrane Flex M	30 mm	20 mm
64052	NeoGen® Collagen Membrane Flex L <i>NeoGen® Collagen Flex resorbs in 3–4 months</i>	40 mm	30 mm
64055	NeoGen® Collagen Membrane Firm S	20 mm	15 mm
64056	NeoGen® Collagen Membrane Firm M	30 mm	20 mm
64057	NeoGen® Collagen Membrane Firm L <i>NeoGen® Collagen Firm resorbs in 6–9 months</i>	40 mm	30 mm

Wound Dressings



Wound Dressings


Art. No.	Description	A	B
64060	NeoGen® Collagen Wound Dressing Tape (1 mm thick, 10 pcs)	75 mm	25 mm
64061	NeoGen® Collagen Wound Dressing Foam (4 mm thick, 10 pcs)	40 mm	20 mm
64062	NeoGen® Collagen Wound Dressing Plug (10 pcs)	20 mm	10 mm



Non-Resorbable Membranes

Ti-reinforced PTFE Membranes

Art. No.	Description	A	B	C	D
64010	NeoGen® PTFE Membrane, Ti-Reinforced – S I (Small Interproximal)	29mm	14mm	19mm	4mm
64011	NeoGen® PTFE Membrane, Ti-Reinforced – M I (Medium Interproximal)	30mm	19mm	19mm	9mm
64012	NeoGen® PTFE Membrane, Ti-Reinforced – L I (Large Interproximal)	36mm	21mm	26mm	11mm
64013	NeoGen® PTFE Membrane, Ti-Reinforced – M (Medium)	32mm	22mm	23mm	12mm
64014	NeoGen® PTFE Membrane, Ti-Reinforced – L (Large)	34mm	25mm	24mm	15mm
64015	NeoGen® PTFE Membrane, Ti-Reinforced – XL (Extra Large)	46mm	32mm	24mm	15mm




Non-reinforced PTFE Membranes

Art. No.	Description	A	B
64019	NeoGen® PTFE Membrane, Non-Reinforced 18×12 mm	18mm	12mm
64020	NeoGen® PTFE Membrane, Non-Reinforced 29×14 mm	29mm	14mm
64021	NeoGen® PTFE Membrane, Non-Reinforced 34×25 mm	34mm	25mm



Fixation Kit

Art. No.	Description
64036	Membrane Tack 3mm – 10 pcs
64031	Tack Positioning Instrument
64033	Tack Mallet
64034	Tack and Screw Cassette
64037	Membrane Screw 3mm – 10 pcs



Sutures


Biotex®

Art. No.	Description
BT3016	USP 3 – 0.16 mm, 3/8 Reverse Cutting Needle 45 cm
BT3019	USP 3 – 0.19 mm, 3/8 Reverse Cutting Needle 45 cm
BT4013	USP 4 – 0.13 mm, 3/8 Reverse Cutting Needle 45 cm
BT4016	USP 4 – 0.16 mm, 3/8 Reverse Cutting Needle 45 cm
BT4019	USP 4 – 0.19 mm, 3/8 Reverse Cutting Needle 45 cm
BT5013	USP 5 – 0.13 mm, 3/8 Reverse Cutting Needle 45 cm
BT5016	USP 5 – 0.16 mm, 3/8 Reverse Cutting Needle 45 cm
BTP401	USP 4 – 0.13 mm, 1/2 Taper Point Needle 45 cm



Other

Art. No.	Description
32.Z4401.00	Neoss PTFE Suture 60 cm 4/0 – Reverse Cutting 16 mm 3/8 (12 pcs)
32.Z4408.00	Neoss PTFE Suture 60 cm 6/0 – Reverse Cutting 12 mm 3/8 (12 pcs)



New opportunities for your digital workflow

Digital libraries can be found under <https://www.neoss.com/en-gb/cad-libraries>



NeoScan

NeoScan

Art. No. Description

91000	NeoScan™ 1000 Intraoral Scanner <i>Includes Handpiece, 4 tips (2×S, 2×L), Holder Set, Manuals, Suitcase. Extended warranties are available for purchase</i>
91012	NeoScan™ 1000 – Tip L TP102N (2 pcs) <i>Includes 2 pcs Tip L</i>
91013	NeoScan™ 1000 – Tip S TP103N (2 pcs) <i>Includes 2 pcs Tip S</i>
91022	Trolley for Intraoral Scanner



Esthetic Healing Abutment with ScanPeg




Esthetic Healing Abutments



ScanPeg



Scan Body



Digital Impression

Esthetic Healing Abutments with ScanPeg

All Esthetic Healing Abutments include Titanium Screw and a ScanPeg

Art. No. Description

		A	B	C
31360	Esthetic Healing Abutment with ScanPeg – Wide Incisor	5mm	5.7mm	6.1mm
31361	Esthetic Healing Abutment with ScanPeg – Narrow Incisor		5.2mm	5mm
31362	Esthetic Healing Abutment with ScanPeg – Canine		6.3mm	5.4mm
31363	Esthetic Healing Abutment with ScanPeg – Pre-molar		6.8mm	5.8mm
31364	Esthetic Healing Abutment with ScanPeg – Molar		7.8mm	7.3mm

ScanPeg

Art. No. Description

Esthetic Healing Abutments with ScanPeg can also be utilised for digital impressions

31365	ScanPeg – 5pcs <i>Used only with Esthetic Healing Abutments</i>
-------	--


Scan Body

Art. No. Description

90181	Scan Body Neoss® SP (IO 1A-A)
90182	Scan Body Neoss® NP (IO 1A-B)
90183	Scan Body Neoss® Access (IO 1B-A)

Analogs

Art. No.	Description
90190	Analog for Printed Model Neoss® SP
90191	Analog for Printed Model Neoss® NP
90192	Analog for Printed Model Neoss® Access
PMA-AIT-1	Elos Accurate® Analog Pliers
PMA-AIP-2	Elos Accurate® Analog Insertion Pin (2 pcs)
PMA-AIS-5	Elos Accurate® Analog Insertion Screw (5 pcs)





Digital Prosthetics

Provisional Abutments

All Provisional Abutments include Provisional Screw and Coping

Art. No.	Description	A	B	C
31267	Provisional Ti Abutment Mono			
31268	Provisional Ti Abutment Multi	15.7 mm	0.7 mm	4.5 mm
31332	Provisional Ti Abutment Mono – Cement Retained			
31241	Provisional Ti Abutment, NP	10.7 mm	0.7 mm	3.8 mm

Includes Provisional Titanium Abutment, Neo Provisional Screw NP and Neo Laboratory Screw NP



NeoBase® Abutments SSC (Straight Screw Channel)

All NeoBase® Abutments SSC include Neo Abutment Screw and Neo Laboratory Screw

Art. No.	Description	A	B	C
31393	NeoBase® Mono SSC G0.3 mm – H3.6 mm	3.9 mm	0.3 mm	
31394	NeoBase® Mono SSC G1.5 mm – H3.6 mm	5.1 mm	1.5 mm	
31395	NeoBase® Mono SSC G0.3 mm – H5.6 mm	5.9 mm	0.3 mm	4.5 mm
31396	NeoBase® Mono SSC G1.5 mm – H5.6 mm	7.1 mm	1.5 mm	
31397	NeoBase® Multi SSC G0.3 mm – H3.6 mm	3.9 mm	0.3 mm	
31398	NeoBase® Mono SSC G0.3 mm – H3.6 mm, NP	3.9 mm		
31399	NeoBase® Mono SSC G0.3 mm – H5.6 mm, NP	5.9 mm	0.3 mm	3.8 mm
31400	NeoBase® Multi SSC G0.3 mm – H3.6 mm, NP	3.9 mm		

Include Neo Abutment Screw and Neo Laboratory Screw





NeoBase® Abutments ASC



NeoBase® Abutments ASC, NP

NeoBase® Abutments ASC (Angulated Screw Channel)

All NeoBase® Abutments ASC include iGO Abutment Screw and iGO Laboratory Screw



Art. No.	Description	A	B	C
31367	NeoBase® Mono ASC G0.3 mm – H3.6 mm	3.9 mm	0.3 mm	4.5 mm
31368	NeoBase® Mono ASC G1.5 mm – H3.6 mm	5.1 mm	1.5 mm	
31369	NeoBase® Mono ASC G0.3 mm – H5.6 mm	5.9 mm	0.3 mm	
31370	NeoBase® Mono ASC G1.5 mm – H5.6 mm	7.1 mm	1.5 mm	
31371	NeoBase® Multi ASC G0.3 mm – H3.6 mm	3.9 mm	0.3 mm	3.8 mm
31372	NeoBase® Mono ASC G0.3 mm – H3.6 mm, NP	3.9 mm		
31373	NeoBase® Mono ASC G0.3 mm – H5.6 mm, NP	5.9 mm	0.3 mm	
31374	NeoBase® Multi ASC G0.3 mm – H3.6 mm, NP	3.9 mm		

Include iGO Abutment Screw NP and iGO Laboratory Screw NP

Neoss TiBase Abutments (compatible with Dentsply Sirona®)

All TiBase Abutments include Neo Abutment Screw and Titanium Screw



TiBase Abutments

Neoss ScanPost

Art. No.	Description	A	B	C
31329	Neoss TiBase N (NB B 3.4 L)	5.3 mm	0.6 mm	4.5 mm
31330	Neoss TiBase W (NB B 4.1 L)			5.0 mm
31345	Neoss TiBase, NP (FX 3.4 S)	5 mm	0.3 mm	3.9 mm
31331	Neoss ScanPost L (TiBase) <i>Includes ScanPost with Screw</i>			
31346	Neoss ScanPost S (TiBase, NP) <i>Includes ScanPost with Screw NP</i>			

Custom Solutions**Neoss NT Blanks****Art. No. Description***Includes Neo Abutment Screw and Neo Laboratory Screw*

31375	Neoss NT Blank Ti SSC, Ø12 mm, SP – 1 pc
31377	Neoss NT Blank Ti SSC, Ø12 mm, NP – 1 pc
31376	Neoss NT Blank CoCr SSC, Ø12 mm, SP – 1 pc
31378	Neoss NT Blank CoCr SSC, Ø12 mm, NP – 1 pc
31391	Neoss NT Blank CoCr SSC, Ø16 mm, SP – 1 pc

Includes Neo Laboratory Screw

31379	Neoss NT Blank Ti SSC, Ø12 mm, SP – 5 pc
31381	Neoss NT Blank Ti SSC, Ø12 mm, NP – 5 pc
31380	Neoss NT Blank CoCr SSC, Ø12 mm, SP – 5 pc
31382	Neoss NT Blank CoCr SSC, Ø12 mm, NP – 5 pc
31392	Neoss NT Blank CoCr SSC, Ø16 mm, SP – 5 pc



NT Blank



ARC

A Neoss Group Company

Neoss Custom Prosthetics – powered by ARC

Order and product information:

Online: www.arcsolutions.se www.neoss-cadcam.de

Phone: +46 42 301 74 40

For full offering, visit www.arcsolutions.se

ARC are specialists in manufacture of custom-made implant prosthetics. With a high focus on quality and staying one step ahead of new technology, ARC places great value in providing customers with maximum service and delivering reliable products in a timely manner.

	Abutment						Bridge				Bar			
	Ti		CoCr		ZrO		Ti		CoCr		Ti		CoCr	
	SSC	ASC	SSC	ASC	SSC	ASC	SSC	ASC	SSC	ASC	SSC	ASC	SSC	ASC
Neoss SP	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Neoss NP	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Neoss Access							✓	✓	✓	✓	✓	✓	✓	✓

SSC= Straight Screw Channel
ASC= Angulated Screw Channel



Full Arch Bridges



Partial Bridges



Custom Bars



Dolder Bars




Vertex Horix Bars



Vertex Bars



Preci classic Bar



SLM Bar and Bridges SSC/ASC



Bar and Bridges ASC



Custom Abutments SSC



Zr solutions on Ti-Base ASC/SSC



Custom Abutments ASC

Design control with *total flexibility*



Esthetic Healing Abutment with ScanPeg

Esthetic Healing Abutments




Healing Abutments Ti

Healing Abutment PEEK with Screw

Wide Healing Abutment

Wide Healing Abutment 5mm



Healing Abutments


Esthetic Healing Abutments with ScanPeg

All Esthetic Healing Abutments include Titanium Screw and a ScanPeg

Art. No.	Description	A	B	C
31360	Esthetic Healing Abutment with ScanPeg – Wide Incisor	5mm	5.7mm	6.1mm
31361	Esthetic Healing Abutment with ScanPeg – Narrow Incisor		5.2mm	5mm
31362	Esthetic Healing Abutment with ScanPeg – Canine		6.3mm	5.4mm
31363	Esthetic Healing Abutment with ScanPeg – Pre-molar		6.8mm	5.8mm
31364	Esthetic Healing Abutment with ScanPeg – Molar		7.8mm	7.3mm

Healing Abutments

Art. No.	Description	A	B
31159	Healing Abutment Ti 2mm	2mm	4mm
31160	Healing Abutment Ti 4mm	4mm	
31161	Healing Abutment Ti 6mm	6mm	
31162	Healing Abutment Ti 8mm	8mm	
31163	Healing Abutment Ti 10mm	10mm	
31148	Healing Abutment PEEK with Screw	5mm	5mm
31320	Healing Abutment PEEK with Screw, NP		4mm
31321	Wide Healing Abutment PEEK with Screw	6.5mm	3.3mm
31335	Wide Healing Abutment PEEK 5mm with Screw	5mm	5mm



Provisional Ti Abutment Mono

Provisional Ti Abutment Multi

Provisional Ti Abutment Mono – Cement Retained

Provisional Ti Abutment, NP

20mm

Coping



Provisional Abutments

Provisional Titanium Abutments

All Provisional Abutments include Provisional Screw and Coping


Art. No.	Description	A	B	C
31267	Provisional Ti Abutment Mono	15.7mm	0.7mm	4.5mm
31268	Provisional Ti Abutment Multi			
31332	Provisional Ti Abutment Mono – Cement Retained	10.7mm	0.7mm	3.8mm
31241	Provisional Ti Abutment, NP			

Includes Provisional Titanium Abutment, Provisional Screw NP and Neo Laboratory Screw NP

Esthetic Tissue Formers

All Esthetic Tissue Formers include Provisional Screw

Art. No.	Description	A	B	C	D
31340	Esthetic Tissue Former – Wide Incisor	10 mm	4 mm	5.7 mm	6.1 mm
31341	Esthetic Tissue Former – Narrow Incisor			5.2 mm	5 mm
31342	Esthetic Tissue Former – Canine			6.3 mm	5.3 mm
31343	Esthetic Tissue Former – Pre-molar			6.8 mm	5.8 mm
31344	Esthetic Tissue Former – Molar			7.8 mm	7.3 mm



Esthetic Tissue Formers




Final Abutments

Prepable Titanium Abutments

All Prepable Abutments include Neo Abutment Screw and Neo Laboratory Screw

Art. No.	Description	A	B	C	D	E
31280	Prepable Ti Abutment 1.5 mm 0° – Wide Incisor	8.5 mm	2.5 mm	1.5 mm	5.8 mm	6.5 mm
31281	Prepable Ti Abutment 1.5 mm 15° – Wide Incisor					
31282	Prepable Ti Abutment 1 mm 0° – Narrow Incisor	8 mm	2 mm	1 mm	5 mm	5.1 mm
31283	Prepable Ti Abutment 1 mm 15° – Narrow Incisor					
31284	Prepable Ti Abutment 3 mm 0° – Narrow Incisor	10 mm	4 mm	3 mm	5.1 mm	5.1 mm
31285	Prepable Ti Abutment 3 mm 15° – Narrow Incisor					
31286	Prepable Ti Abutment 3 mm 0° – Canine	10 mm	4 mm	3 mm	6.7 mm	5.4 mm
31305	Prepable Ti Abutment 1.5 mm 0° – Pre-molar	7 mm	2 mm	1.5 mm	6.9 mm	5.7 mm
31287	Prepable Ti Abutment 3 mm 0° – Pre-molar	9 mm	3 mm	3 mm	6.9 mm	5.7 mm
31306	Prepable Ti Abutment 1.5 mm 0° – Molar	7 mm	2 mm	1.5 mm	7.8 mm	7.1 mm
31288	Prepable Ti Abutment 3 mm 0° – Molar	9 mm	3 mm	3 mm	7.8 mm	7.1 mm
31201	Prepable Ti Abutment 1 mm 0°, NP	8 mm	2 mm	1 mm	4 mm	
31202	Prepable Ti Abutment 1 mm 15°, NP					



Wide Incisors

Narrow Incisors



Canine, Pre-molar and Molar

Prepable Ti Abutments

Prepable Titanium Abutments, Alternative Emergence Profiles

Art. No.	Description	A	B	C	D	E
31317	Prepable Ti Abutment 2 mm 0° – Wide Emergence <i>Includes Neo Abutment Screw, Neo Laboratory Screw, specific Cover Screw and specific Healing Abutment PEEK with Screw. Compatible with Implant diameters Ø5.0 – 6.5 mm</i>	9 mm	2 mm	2 mm	7.9 mm	7.1 mm
31318	Prepable Ti Abutment 1.5 mm 0° – Narrow Emergence	8.5 mm	2 mm	1.5 mm	4.4 mm	4.5 mm
31319	Prepable Ti Abutment 1.5 mm 15° – Narrow Emergence <i>All Narrow Emergence Prepable Abutments include Neo Abutment Screw and Neo Laboratory Screw. Only use with Ø3.5 and Ø4.0 mm implants</i>					



Wide Emergence



Narrow Emergence




Zirconia Abutments

Zirconia Abutments

All Zirconia Abutments include separate Zirconia Coping, Ti NeoLink® Mono, Neo Abutment Screw and Neo Laboratory Screw

Art. No.	Description	A	B	C	D	E
31290	Zirconia Abutment – Wide Incisor	9mm	2.5mm	1.5mm	5.8mm	6.4mm
31291	Zirconia Abutment – Narrow Incisor				5.4mm	5.1mm
31292	Zirconia Abutment – Canine		2mm	2mm	6.2mm	5.8mm
31293	Zirconia Abutment – Pre-molar				6.4mm	5.5mm
31294	Zirconia Abutment – Molar				7.2mm	6.7mm




Final Restorations – Casting

Gold NeoLink® (Direct Casting)

Gold NeoLinks include Neo Abutment Screw, Neo Laboratory Screw and Straight Copings



Gold NeoLink®



19mm

Copings

Art. No.	Description	A	B	C
31102	Gold NeoLink® Mono	4.0mm	0.4mm	4.5mm
31103	Gold NeoLink® Multi			
31199	Gold NeoLink® Mono, NP <i>Includes Neo Abutment Screw Ti NP, Neo Laboratory Screw NP and a set of Preformed Plastic Copings</i>			3.8mm

CoCr Abutments (Indirect Casting)

CoCr Abutments include Neo Abutment Screw, Neo Laboratory Screw and Straight Coping



CoCr Abutment Mono



CoCr Abutment Multi



Burnout Abutment Mono



Burnout Abutment Multi



Access Burnout Abutment Multi



Art. No.	Description	A	B	C
31269	CoCr Abutment Mono	15.7mm	0.7mm	4.5mm
31260	CoCr Abutment Multi			

Burnout Abutments (Direct Casting)

31333, 31157 and 31266 include Straight Coping


Art. No.	Description	A	B	C
31333	Burnout Abutment Mono – 2 pcs	15.7mm	0.7mm	4.5mm
31157	Burnout Abutment Multi – 2 pcs			
31266	Access Burnout Abutment Multi – 2 pcs			
31255	Burnout NeoLink® Multi NP – 2 pcs <i>Includes set of Preformed Plastic Copings</i>	4.3mm		3.8mm



Burnout NeoLink® Multi



NeoLink® Plastic Copings



Ti NeoLink®


Titanium NeoLinks include Neo Abutment Screw, Neo Laboratory Screw and Straight Copings

Art. No.	Description	A	B	C
31133	Ti NeoLink® Mono			4.5mm
31134	Ti NeoLink® Multi	3.9mm	0.3mm	
31234	Ti NeoLink® Mono, NP			3.8mm

Includes Neo Abutment Screw Ti NP, Neo Laboratory Screw NP and a set of Preformed Plastic Copings



Ti NeoLink®



Multiple Unit Abutments

Access Abutments (Screw-retained multiple unit restoration)

Art. No.	Description	A	B	C
<i>All Access Abutments include Abutment and a Healing Abutment</i>				
31238	Access Abutment 1.5 mm	1.5mm	4 mm	
31215	Access Abutment 3.0 mm	3 mm		
31216	Access Abutment 4.0 mm	4 mm		
31218	Access Abutment 6.0 mm	6 mm		
31219	Access Abutment 10° 2.6 mm	2.6mm	4 mm	
31243	Access Abutment 10° 4.6 mm	4.6mm		
31221	Access Abutment 20° 2.6 mm	2.6mm		
31244	Access Abutment 20° 4.6 mm	4.6mm		
31223	Access Abutment 30° 2.9 mm	2.9mm		
31245	Access Abutment 30° 4.9 mm	4.9mm		
31225	Access Impression Coping 8 mm and Replica <i>Includes Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>			
31366	Access Impression Coping 13 mm and Replica			
90183	Scan Body Neoss® Access (IO 1B-A)			
31226	Access Ti NeoLink® Multi	4.3mm	0.3mm	4.5mm
31227	Access Gold NeoLink® Multi <i>31226 and 31227 include NeoLink®, Prosthetic Screw, Provisional/Laboratory Screw and a Plastic Coping</i>	4.7mm	0.7mm	4.5mm



Access Abutments



Access Impression Coping and Replica


Scan Body



Access Ti NeoLink®



Access Gold NeoLink®



Art. No.	Description	A	B	C
31334	Access Ball Abutment 1mm			
90207	Access Equator Abutment Kit 0.5 mm			
90208	Access Equator Abutment Kit 1.5 mm			
31301	Access Provisional Ti Abutment Multi <i>Includes Provisional/Laboratory Screw and Straight Coping</i>	15.7 mm	0.7 mm	4.5 mm
31266	Access Burnout Abutment Multi – 2 pcs <i>Includes Straight Coping</i>			
31229	Access Neo Prosthetic Screw – 1 pc			
31257	Access Neo Provisional/Laboratory Screw – 5 pcs			
31265	Access Neo Laboratory Screw Long – 5 pcs			
31231	Access Replica – 5 pcs			
90192	Analog for Printed Model Neoss® Access			

Overdenture Abutments

Equator Abutments

All kits include Equator Abutment, Housing with processing Cap, Protector Disk, Standard, Strong, Soft and Extra Soft Retention caps




Art. No.	Description
90201	Neoss Equator Abutment Kit 1mm
90202	Neoss Equator Abutment Kit 2mm
90242	Neoss Equator Abutment Kit 2mm, NP
90203	Neoss Equator Abutment Kit 3mm
90243	Neoss Equator Abutment Kit 3mm, NP
90204	Neoss Equator Abutment Kit 4mm
90244	Neoss Equator Abutment Kit 4mm, NP
90205	Neoss Equator Abutment Kit 5mm
90206	Neoss Equator Abutment Kit 6mm
90207	Access Equator Abutment Kit 0.5mm
90208	Access Equator Abutment Kit 1.5mm
90210	Neoss Equator Cap and Housing Kit
90215	Neoss Equator Cap and Housing Kit – Extended Divergence
90211	Neoss Equator Cap Standard Retention, White – 2 pcs
90212	Neoss Equator Cap Strong Retention, Violet – 2 pcs
90213	Neoss Equator Cap Soft Retention, Pink – 2 pcs
90214	Neoss Equator Cap Extra Soft Retention, Yellow – 2 pcs
90220	Neoss Equator Impression Coping & Replica – 2 sets
90222	Neoss Equator Protector Disk – 10 pcs
90232	Neoss Equator Cap Tool

Ball Abutments


All Ball Abutments include Replica

Art. No.	Description
31140	Ball Abutment 1mm
31141	Ball Abutment 3mm
31142	Ball Abutment 4.5mm
31143	Ball Abutment 6mm
31334	Access Ball Abutment 1mm
51124	Ball Driver
SF-434-9-A	LV SFE Female yellow – 6 pcs
SF-434-9-B	LV SFE Female white – 6 pcs
SF-434-9-C	LV SFE Female red – 6 pcs
AC-474-1	LV Insertion Tool/0.9 mm hex driver
SF-474-Q	LV SFE Space Maintenance – 6 pcs
SF-444-3TI	LV SFE Housing Ti – 6 pcs



Bar Abutment


Art. No.	Description
31168	Bar Abutment Gold <i>Includes Abutment, Neo Abutment Screw and Neo Laboratory Screw</i>
31261	Bar Abutment Ti 14 mm <i>Includes Neoss Abutment Screw and Laboratory Screw</i>



Impressions & Replicas


Impressions Copings

Art. No.	Description
31170	Impression Coping 8mm and Replica
31145	Impression Coping 11mm and Replica
31183	Impression Coping 18mm and Replica <i>31170, 31145 and 31183 include Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>
31302	Impression Coping Open Tray 8mm and Replica
31185	Impression Coping Open Tray 13mm and Replica <i>31302 and 31185 include Coping, Screw, Replica and Extension Tube</i>
31197	Impression Coping 11mm and Replica, NP <i>Includes Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>
31366	Access Impression Coping 13mm and Replica
31186	Impression Coping Cap Closed Tray – 10 pcs
31262	Impression Coping Extension Tube – 20 pcs



Replicas

Art. No.	Description
31231	Access Replica – 5pcs
31146	Implant Replica – 5pcs
31198	Implant Replica, NP – 5pcs
31304	Implant Replica, NP – 1pc
31166	Protection Replica – 1pc





Screws

Art. No.	Description
31258	Neo Abutment Screw, SP – 1pc
31259	Neo Abutment Screw, SP – 5pcs
31240	Neo Abutment Screw Ti, NP – 1pc
31149	Neo Laboratory Screw, SP – 5pcs
31239	Neo Laboratory Screw, NP – 1pc
31303	Neo Laboratory Screw Long, SP – 5pcs
31347	iGO Abutment Screw, SP – 1pc
31349	iGO Abutment Screw, NP – 1pc
31350	iGO Laboratory Screw, SP – 1pc
31352	iGO Laboratory Screw, NP – 1pc
31229	Access Neo Prosthetic Screw – 1pc
31257	Access Neo Provisional/Laboratory Screw – 5pcs
31265	Access Neo Laboratory Screw Long – 5pcs

Neo = Neoss standard
Screw / Screwdriver



The instruments *you need*

Drills

Art. No. Description

41153	Round Bur Ø1.8 mm
41152	Lance Drill Ø1.5×3/5/7 mm
41170	Pilot Drill Ø2.2/3.0 mm
41154	Twist Drill, Straight Ø2.2×7 – 17 mm
41156	Twist Drill, Straight Ø2.85×7 – 17 mm
41157	Twist Drill, Straight Ø3.0×7 – 17 mm
41158	Twist Drill, Straight Ø3.2×7 – 17 mm
41159	Twist Drill, Straight Ø3.4×7 – 17 mm
41160	Twist Drill, Straight Ø3.6×7 – 17 mm
41161	Twist Drill, Straight Ø3.9×7 – 17 mm
41162	Twist Drill, Straight Ø4.1×7 – 17 mm
41163	Twist Drill, Straight Ø4.4×7 – 17 mm
41164	Twist Drill, Straight Ø4.6×7 – 17 mm
41165	Twist Drill, Straight Ø4.9×7 – 17 mm
41166	Twist Drill, Straight Ø5.1×7 – 17 mm

41205 Twist Drill, Straight Ø6.0×7 – 11 (13) mm

41193 Drill Kit, Straight Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2, 3.0, 3.4, 3.9 and 4.4

41172	Twist Drill, Tapered Ø3.0×7 – 17 mm	1.9 mm
41173	Twist Drill, Tapered Ø3.4×7 – 17 mm	2.3 mm
41174	Twist Drill, Tapered Ø3.9×7 – 17 mm	2.8 mm
41175	Twist Drill, Tapered Ø4.4×7 – 17 mm	3.3 mm
41176	Twist Drill, Tapered Ø4.9×7 – 17 mm	3.8 mm
41206	Twist Drill, Tapered Ø5.5×7 – 17 mm	4.3 mm
41204	Twist Drill, Tapered Ø5.5×7 – 11 (13) mm	4.3 mm

41192 Drill Kit, Tapered Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2 S, 3.0 T, 3.4 T, 3.9 T and 4.4 T

41177 Drill Kit, Edge Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2 S, 3.0 T, 3.4 T, 3.9 T and 4.4 T

41142	Drill Stop 7 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41143	Drill Stop 9 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41144	Drill Stop 11 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41145	Drill Stop 13 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41146	Drill Stop 15 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm

41215	Twist Drill, Straight Ø2.2×7 – 17 mm
41216	Twist Drill, Straight Ø2.85×7 – 17 mm
41217	Twist Drill, Straight Ø3.0×7 – 17 mm
41218	Twist Drill, Straight Ø3.2×7 – 17 mm
41219	Twist Drill, Straight Ø3.4×7 – 17 mm
41220	Twist Drill, Straight Ø3.6×7 – 17 mm
41221	Twist Drill, Straight Ø3.9×7 – 17 mm
41222	Twist Drill, Straight Ø4.1×7 – 17 mm
41223	Twist Drill, Straight Ø4.4×7 – 17 mm
41224	Twist Drill, Straight Ø4.6×7 – 17 mm
41225	Twist Drill, Straight Ø4.9×7 – 17 mm
41226	Twist Drill, Straight Ø5.1×7 – 17 mm



Round Bur



Lance Drill



Pilot Drill



Twist Drill, Straight



Twist Drill, Tapered



Drill Stop

Single Use

Multiple Use



Twist Drill, Straight – Short



Twist Drill, Tapered – Short



Countersink



Countersink, Straight




Countersink, Tapered



Countersink, Edge



Screw Tap



Bone Profiler



Bone Remover Cover Screw

Art. No. Description

41233 Drill Kit, Straight Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2, 3.0, 3.4, 3.9 and 4.4

41184 Drill Kit, Straight – Short 7 – 13 mm
Includes Twist Drills Ø2.2, 3.0, 3.4, 3.9, 4.4 and 4.9

41227 Twist Drill, Tapered Ø3.0×7 – 17 mm

41228 Twist Drill, Tapered Ø3.4×7 – 17 mm

41229 Twist Drill, Tapered Ø3.9×7 – 17 mm

41230 Twist Drill, Tapered Ø4.4×7 – 17 mm

41231 Twist Drill, Tapered Ø4.9×7 – 17 mm

41232 Twist Drill, Tapered Ø5.5×7 – 17 mm

41235 Drill Kit, Tapered Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T and 4.4T

41191 Drill Kit Tapered – Short 7 – 13 mm
Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T, 4.4T, and 4.9T

41236 Drill Kit, Edge Implants Ø3.5 – 5.0
Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T and 4.4T

41139 Countersink, Straight Ø3.25

41110 Countersink, Straight Ø3.5

41111 Countersink, Straight Ø4.0

41112 Countersink, Straight Ø4.5

41133 Countersink, Straight Ø5.0

41127 Countersink, Straight Ø5.5

41178 Countersink, Tapered Ø3.5

41179 Countersink, Tapered Ø4.0

41180 Countersink, Tapered Ø4.5

41181 Countersink, Tapered Ø5.0

41182 Countersink, Tapered Ø5.5

41210 Countersink, Edge Ø3.5

41211 Countersink, Edge Ø4.0

41212 Countersink, Edge Ø4.5

41213 Countersink, Edge Ø5.0

41202 Countersink Ø6.0

41201 Countersink Ø6.5

41113 Screw Tap Ø3.5

41114 Screw Tap Ø4.0

41115 Screw Tap Ø4.5

41134 Screw Tap Ø5.0

41126 Screw Tap Ø5.5

41140 Screw Tap Ø3.25

41237 Bone Profiler Ø5.1mm, SP

41238 Bone Profiler Ø6.6mm, SP

41239 Bone Profiler Ø4.1 mm, NP

41240 Bone Remover Cover Screw, SP

41241 Bone Remover Cover Screw, NP

19mm
2.3mm
2.8mm
3.3mm
3.8mm
4.3mm

Ø
2.7mm
2.9mm
3.3mm
3.5mm
4.3mm
4.7mm
2.0mm
2.8mm
3.3mm
3.7mm
4.2mm
1.8mm
2.6mm
3.0mm
3.5mm
4.8mm
5.2mm

Multiple Use

Kits


Surgical and Prosthetic Trays

Art. No. Description

51185 Neoss Clinical Organiser, Edge

51150 Neoss Clinical Organiser, Straight


51151 Neoss Clinical Organiser, Tapered



Clinical Organiser, Edge



Clinical Organiser, Tapered



Clinical Organiser, Straight

Trays & Instrument Kits

Art. No. Description

51180	Neoss Prosthetic Tray – Neo
51181	Neoss Prosthetic Instrument Kit – Neo <i>Includes Prosthetic Tray, Ratchet, Screwdriver Machine 22mm, Screwdriver Machine 32mm, Screwdriver Manual 15mm and Manual Handle</i>
41151	Surgical Instrument Kit <i>Includes art no 41120, 51118, 51119, 51144, 51139, 51140, 51125, 51126</i>
51131	Neoss System Tray <i>Includes your choice of 51150, 51151 or 51185</i>



Neoss Prosthetic Tray – Neo



Neoss Prosthetic Instrument Kit – Neo



Neoss System Tray



Neoss System Tray – Straight

Guide

Art. No. Description

51154	Neoss® Guide Kit Ø3.5 – Ø6.0 <i>(Neoss Implants Ø3.5 – Ø6.0mm, 7 – 15mm lengths.) Includes 5 × Drill-Hubs, 2 × Countersink Conduits, 4 × Implant Mount Ø5.0 S9, 2 × Implant Mount Ø6.0 S9, 7 × Guide Keys, 3 × Guide Key Adapter, 1 × Holder Guide Key Adapter. Spare parts available, contact local representative for details</i>
51173	Neoss® Guide – Implant Mount SP, Ø5.0 S9 – 2 pcs
51174	Neoss® Guide – Implant Mount SP, Ø6.0 S9 – 2 pcs
51156	Neoss® Guide – Sleeve Ø5.0 – 10 pcs
51157	Neoss® Guide – Sleeve Ø6.0 – 5 pcs
51155	Neoss® Guide – Implant Inserters, SP <i>For partially guided only, see IFU</i>



Neoss® Guide Kit




Implant Mount



Sleeve



Implant Inserters



Drill-Hubs



Guide Key



Countersink Conduit



Drill Sleeve



Guide Key Adapter



Holder Guide Key Adapter

Neoss® Guide Kit – SPARE PARTS

Art. No. Description

51160	Neoss® Guide – Drill-Hub I
51161	Neoss® Guide – Drill-Hub II
51162	Neoss® Guide – Drill-Hub III
51163	Neoss® Guide – Drill-Hub IV
51164	Neoss® Guide – Drill-Hub V
51165	Neoss® Guide – Guide Key Ø2.2/2.85
51166	Neoss® Guide – Guide Key Ø3.0/3.2
51167	Neoss® Guide – Guide Key Ø3.4/3.6
51168	Neoss® Guide – Guide Key Ø3.9/4.1
51169	Neoss® Guide – Guide Key Ø4.4/4.6
51170	Neoss® Guide – Guide Key Ø4.9/5.1
51178	Neoss® Guide – Guide Key Ø5.5
51171	Neoss® Guide – Countersink Conduit Ø5.0
51172	Neoss® Guide – Countersink Conduit Ø6.0
51176	Neoss® Guide – Guide Key Adapter – 2 pcs
51177	Neoss® Guide – Holder Guide Key Adapter
51158	Drill Sleeve Ø2.2 – 10 pcs


Demo Kit

Art. No. Description

61150	King Size Dummy Kit <i>Includes Neoss Straight, Tapered and Edge Implants, implant inserter, screwdriver, abutment screw, NeoLink® Abutment, Esthetic Abutment Narrow Incisor with ScanPeg and crown</i>
-------	---



King Size Dummy Kit



Instruments

Art. No.	Description
51137	Implant Inserter, SP 17 mm
51118	Implant Inserter, SP 22 mm
51119	Implant Inserter, SP 32 mm
51145	Implant Inserter, NP 24 mm
51142	Implant Inserter, NP 32 mm
51139	Neo Screwdriver, Machine 22 mm
51140	Neo Screwdriver, Machine 32 mm
51141	Neo Screwdriver, Manual 15 mm
51183	iGO Screwdriver, Machine 22 mm
51184	iGO Screwdriver, Machine 32 mm
51125	Direction/Depth Gauge – 4 pcs
41120	Drill Extender
51121	Ratchet
51126	Manual Handle
51144	Wrench Adapter
51149	Angulation Gauge, SP
90138	Adjustment Handle
41214	X-ray Planner ProActive Edge Implants
41119	X-ray Planner ProActive Straight Implants
41183	X-ray Planner ProActive Tapered Implants
41203	X-ray Planner Ø6.0 Implants
41200	X-ray Planner Ø6.5 Implants

Neo = Neoss standard
Screw/Screwdriver



Surgical Equipment

Neoss Penguin^{RFA}

Art. No.	Description
90300	Neoss Penguin ^{RFA} Instrument <i>Includes Neoss Penguin^{RFA} Instrument, Neoss MulTipegTM Driver, charger and adapters</i>
90307	Sterile Cover Neoss Penguin ^{RFA} Instrument – 20 pcs
90306	Neoss MulTipeg TM Driver
90310	MulTipeg TM 18, Neoss SP
90311	MulTipeg TM 19, Neoss NP
90312	MulTipeg TM 52, Neoss Access

MulTipegTMs for other brands are available, contact your local representative for details



Prosthetic Solution Standard and Narrow Platform, overview

Overview Prosthetic Solutions, SP and NP (Ø3.25 – 6.5 Implants)



		Units		Retention		Implant platform		Level		Material					Angulated screw channel	Esthetiline
		Mono	Multi	Cement retained	Screw retained	SP	NP	Implant level	Access level	Ti	ZrO	Gold	CoCr	Polymer		
Provisional Restorations	Provisional Ti Abutment	✓	✓	✓	✓	✓	✓	✓	✓	✓						
	Esthetic Tissue Former	✓		✓	✓	✓		✓						✓		✓
Final Restorations	Prefabricated	Prepable Ti Abutment	✓		✓	✓	✓	✓		✓						✓
		Ti NeoLink®	✓	✓	✓	✓	✓	✓	✓	✓						✓
		Zirconia Abutment	✓		✓	✓		✓			✓					✓
		Access Abutment (0 – 30°)		✓	✓	✓		✓		✓						
	Digital	NeoBase® (SSC, ASC*)	✓	✓	✓	✓	✓	✓		✓					✓	✓
		TiBase	✓		✓	✓	✓	✓		✓						
		Blanks	✓		✓	✓	✓	✓		✓			✓			✓
		Individual Prosthetics* central milling	✓	✓	✓	✓	✓	✓	✓	✓	✓		✓		✓	✓
	Casting	Gold NeoLink®	✓	✓		✓	✓	✓	✓			✓				✓
		Burnout Abutment	✓	✓	✓	✓	✓	✓	✓							
		CoCr Abutment	✓	✓	✓	✓		✓					✓			
Overdenture		Equator Abutment		✓		✓	✓	✓	✓							
		Ball Abutment		✓		✓		✓								
		Bar Abutment		✓		✓		✓								

* Utilise iGO Abutment Screw and Screwdriver for angulated screw channels, optional for straight screw channels on Individual Prosthetics milled by Neoss

Torque and Speed recommendation, overview

Neoss Implant System Torque Recommendation (Ncm)

Implants	Screw Taps	Healing			Provisional		Impression	Final restoration					
		Cover Screw	Healing Abutment PEEK and Esthetic Healing Abutment	Healing Abutment Ti	Esthetic Tissue Former	Provisional Ti Abutment	Impression Coping Screw	Neo Abutment Screw	Access Abutment Straight & Angulated	Neo Access Prosthetic Screw	Ball Abutment including Access level	Equator Abutment	
45 Max	40 – 45	10 Max	10 Max	10 Max	20	20	10 Max	32*	32	20	20	Implant level	Access level
												32	20

* 35 Ncm optional for high-load cases

Final restoration with iGO Screws		
iGO Abutment Screw, SP	iGO Abutment Screw, NP	iGO Prosthetic Screw, Access
32	20	20

Neoss Implant System Drilling/Insertion Speed Recommendation (rpm)

Drills	Screw Taps	Countersinks	Implants	Bone Mill
800 – 2000**	20	800	20 Max	40 Max

** lower speed for larger drills

Symbols, names and abbreviations



SP = Standard Platform Ø3.5 – 6.5 Implants



NP = Narrow Platform Ø3.25 Implants



iGO = utilises iGO Screw and Screwdriver

Neo = Neoss standard Screw/Screwdriver

SSC = Straight Screw Channel

ASC = Angulated Screw Channel

You can now buy your products online!

Products can also be ordered online at **neoss.com**.

Some products may not be regulatory or commercially available for all markets and available in limited quantities.

Contact your local representative for availability and prices.

Disclaimer of Liability

Neoss products may only be used according to the manufacturers' instructions and recommendations.

The user of Neoss products should determine their suitability for particular patients and indications.

Neoss Limited disclaims any liability, expressed or implied, and shall have no responsibility for any direct, indirect, punitive or other damages arising out of or in connection with any errors in professional judgment or practice in the use or placement of the Neoss products.

Product images are not to scale and are subject to change.

All dimensions are for identification purposes only.

Copyrights, design rights and trademarks:

Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorisation of Neoss Limited.

Neoss, 'N' logo, NeoGen, NeoLink, NeoLoc, NeoBase and Neoss ProActive are registered trademarks of Neoss Group.

NeoScan is a trademark of Neoss Group.

THE Graft™ is a trademark of Purgo Biologics Inc. Biotex® is a registered trademark by Purgo Biologics Inc.

Elos Accurate is a registered trademark of Elos Medtech AB.



Contact details

Our offices

ARC SOLUTIONS

Arc Solutions
Box 13011
SE-250 13 Helsingborg
Sweden
T: +46 42 301 74 40
E: info@arcsolutions.se

AUSTRALIA & NEW ZEALAND

Neoss Australia Pty. Ltd
PO Box 404
New Farm
QLD 4005
Phone: +61 7 3216 0165
E: info.au@neoss.com

CHINA

Neoss China
7F, #303, Songhu Rd. (KIC Build #11)
Yangpu District
Shanghai, 200433
PR China
T: +86 21 80143135
E: china@neoss.com

EUROPEAN DISTRIBUTION

Neoss AG
The Circle 23
8058 Zurich-Airport
Switzerland
E: distributors@neoss.com

GERMANY & AUSTRIA

Neoss GmbH
Im Mediapark 5b
50670 Köln
T: +49 221 96980 10
F: +49 221 96980 199
E: info@neoss.de

ITALY

Neoss Italia S.r.l.
Viale Certosa 138
IT-20156 Milano
T: +39 02 92952 1
E: italia@neoss.info

MIDDLE EAST & AFRICA

Neoss Limited (DMCC BRANCH)
Unit 1505 Fortune Tower
Cluster C – JLT
Dubai
P.O. Box 64093
United Arab Emirates
E: info.mea@neoss.com

SWEDEN, DENMARK & NORWAY

Neoss AB
Arvid Wallgrens backe 20
SE-413 46 Göteborg
T: +46 31 88 12 80
E: info@neoss.se

UNITED KINGDOM & REPUBLIC OF IRELAND

Neoss Ltd
Windsor House
Cornwall Road
Harrogate
HG1 2PW
T: +44 1423 817 733
F: +44 1423 817 744
E: info@neoss.com

UNITED STATES

Neoss Inc.
21860 Burbank Blvd.,
Ste. 190 (North Entrance)
Woodland Hills, CA 91367
T: +1 866 626 3677
F: +1 818 432 2640
E: contact.usa@neoss.com

Our distributors

BAHRAIN

Cigalah Gulf Medical
Building 670, Block 309, Road 919
Salmaniya, Bahrain
T: +973 17232170
E: support@cgmed.com
maha@cgmed.com

CROATIA

Medical Intertrade d.o.o
10431 Sveta Nedelja
HRVATSKA
T: +385 1 33 74 020
F: +385 1 33 36 035
E: dental@medical-intertrade.hr

EGYPT

Dental Future
Building No 5
Mohamed Hassanein Heikal
St. Off Abbas El Akkad Street
1st Floor Flat 3
Cairo, Egypt
T: +20 109 6696449
E: mahmod@dentalfuturekw.com

IRAQ

Eblal Scientific Bureau
Al-Ameerat street
Almansour
Baghdad
Iraq
T: +964 780 411 0110
E: info@eblal.com

JAPAN

Neoss Japan Co.
Y101 Building 5F
No. 1-4 Kanda-Sudacho Chiyodaku
Tokyo 101-0041
T: +81 3 5289 3511
F: +81 3 5289 3512
E: info@neoss.co.jp

JORDAN

Advanced Intelligence for Medical Supplies & Equipment Ltd.
Al-Swefeyyeh Paris str
Alhusainy Complex
Amman 11196
Jordan
T: +962 7 97389535
F: +962 6 5826797
E: f.khaddam@advanceime.com

KINGDOM OF SAUDI ARABIA

Tamer Group
Mohammad Ibn Abdul Wahab Street
Al Ammariya District
P.O. Box 180
Jeddah, 21411
Kingdom of Saudi Arabia
T: +966 12 643 5600
F: +966 12 642 4141
E: alaa.wanly@tamergroup.com

KUWAIT

Dental Future Company
Salimiya Salem el Mubark street
Suhaib complex office 10
Kuwait City
Kuwait
T: +965 51155997
E: Info@dentalfuturekw.com

LEBANON

Dentaltech Sarl
Brasilia Street
Hikmeh Building
First Basement
Baabda, Beirut – Lebanon
T: +961 5 950 707
F: +961 5 950 808

NETHERLANDS

Rotor, Medical Supplies B.V.
Rijksweg Zuid 51
NL-6071HV Swalmen
T: +31 46 481 12 02
F: +31 46 481 12 04
E: mail@rotor-medical.com

OMAN

Al Zahrawi Medical Supplies LLC.
Office No: 16 , Al Dar building opposite
Al Zubair automotive, street 61, Ghala
area.
Muscat, Oman
P.O. Box: 1627
T: +968 24587308
F: +968 24587487
E: info.oman@zahrawigroup.com

QATAR

Al Zahrawi Medical LLC
Building No. 36
Zone No. 39
Al Sadd Area
Street No.808
Doha, Qatar, P.O. Box 63065
T: +974 4431 872
F: +974 44365216
E: info.qatar@zahrawigroup.com

SPAIN

Formación, Gestión Y Rendimiento, S.L.
C. Nil Fabra 34-36.
Etslo 2 B
08012, Barcelona
Spain
T: +34 931 408 770
E: info@fgdental.com

SWITZERLAND

Flexident AG
Amlehnstr. 22
CH-6010 Kriens
T: +41 41 3104020
E: contact@flexident.ch

THAILAND

S.D. Tantavade Co., Ltd.
47 Samakkee 20, Samakkee Road
Thasai Muang
Nonthaburi 11000
Thailand
T: +66 2 952 4815 18
F: +66 2 952 4819
E: info@sdt1988.com

TURKEY

Dentmark Sağlık Ürünleri Anonim Şirketi
Cevizli Mah. Cevizli D-100 Güney
Yanyol Cad.
No:2 Kuruş Kule Kat:3 D:52
Kartal-İstanbul PO: 34865
T: +90 216 417 4040
E: info@dentmark.com.tr

UNITED ARAB EMIRATES (ABU DHABI)

Elwan Technical Supplies
Airport Street Abu Dhabi UAE
T: +971 26393292
E: info@elwandental.com

UNITED ARAB EMIRATES (DUBAI & NORTHERN EMIRATES)

Al Zahrawi Medical Supplies
1st Floor, Ford Middle East & Africa Building
Al Wasl Office Complex, Sheikh Zayed Road, Downtown
Dubai, United Arab Emirates
P.O. Box: 5973
T: +971 4 262 2728
F: +971 4 262 5506
E: info@zahrawigroup.com

For our Terms & Conditions consult your local website on www.neoss.com

Contact your local representative for product availability details.
All products are not available in every market.