



Off-the-shelf revolutionary *2-1 solution*

Esthetic Healing Abutments with ScanPeg

Intelligent Simplicity

A smarter, *more efficient patient-friendly* process

The Neoss® Esthetic Healing Abutments with ScanPeg is a combination of healing abutment and scan abutment in a revolutionary 2-in-1 solution. This is a patient-friendly, simplified and effortless procedure for accurate single tooth intraoral scanning and reduces the risk of soft-tissue damage with multiple healing abutment removals.

- ✓ Patient friendly
- ✓ Less invasive
- ✓ Easy handling
- ✓ Excellent esthetics



Providing treatment with *fewer steps*

Neoss always strives at developing products that provide a treatment that is less invasive and can increase patient comfort. The Neoss Esthetic Healing Abutments with ScanPeg is part of the NeossONE solution, a prosthetic platform where the same abutment fits all Neoss implants no matter the implant design or diameter. Providing a complete range of anatomically shaped PEEK Healing Abutments with a dedicated Scan Abutment, for simplified and accurate intraoral scanning.

Easier handling provides, faster improved procedure

- Minimized component exchange
- Streamlined intraoral scanning
- Less chair time, fewer treatment steps



More patient friendly due to undisturbed healing

- Biological seal is maintained, tissue level is preserved
- Less invasive for increased patient comfort

Providing excellent esthetics

- Shapes soft tissue during healing
- Perfectly matched with Neoss Esthetiline final abutments and Neoss customized abutments

A step ahead

01 Healing abutment placement


The Esthetic Healing Abutment is placed during first or second stage surgery and then left until the final abutment is placed.



02 ScanPeg placement and scanning

An accurate intraoral scan is carried out at healing abutment level utilizing the ScanPeg. The ScanPeg comes co-packed with the Esthetic Healing Abutment and is easy to both attach and remove.

The intraoral digital impression can be taken whenever you choose. ScanPeg shape is included in Neoss CAD libraries.



03 Placement of the final restoration

Flexibility in when you take the digital impression allows design work to continue during the healing phase. The healing abutment is not removed until the final restoration is placed, eliminating one complete surgical intervention.


Contours of Esthetic Healing Abutments are included in Neoss CAD libraries to design the emergence profile from patients contour to final restoration.



exocad 3shape

Assortment

Article no.	Description
31360	Esthetic Healing Abutment with ScanPeg – Wide Incisor
31361	Esthetic Healing Abutment with ScanPeg – Narrow Incisor
31362	Esthetic Healing Abutment with ScanPeg – Canine
31363	Esthetic Healing Abutment with ScanPeg – Premolar
31364	Esthetic Healing Abutment with ScanPeg – Molar
31365	ScanPeg – 5 pcs



Part of the Neoss Esthetiline range

Hear it from the *experts*



“With the ScanPeg, a surgeon or prosthodontist can make a digital impression at any point in time without disturbing the soft-tissue healing. *It is more than simple, and absolutely intelligent.*”

– Dr. Jakob Zwaan
Dentist | Calusco D'Adda, Italy



“The NeoScan 1000 is *very fast, easy to use, and delivers extremely precise data.* Together with the Esthetic Healing Abutment with ScanPeg, the digital workflow is raised to a new level: *Greater patient comfort, more time efficiency and undisturbed healing for better esthetics and function.*”

– Dr. Fabia Flensburg
Dentist | Essen, Germany



“The Esthetic Healing Abutment has completely transformed the way in which we can control the emergence healing of a recently placed implant with *minimal effort and maximum results.* This is the only fully integrated system made specifically for digital dentistry, the ScanPeg truly allows for *one abutment, one time concept.* This is the best implementation of a restorative system in dentistry right now.”

– Dr. Robert Ritter
Dentist | Florida, United States



“The Esthetic Healing Abutments with ScanPeg are an excellent addition to the Neoss range. *They simplify the restorative process* and ensure you don't need to disrupt the soft tissue. The ability of the lab to copy the esthetic emergence profile makes fitting restorations simple and *reduces excessive pressure on the soft tissues.*”

– Dr. Deepesh Patel
Specialist in prosthodontics | Cambridge, United Kingdom

Your solution for *digital impressions*

The NeoScan™ 1000 is a fast, precise and easy-to-use intraoral scanner that provides the possibility for a flexible workflow. The simplified USB cable connection gives you the option to share the scanner between dental chairs in your practice. No separate power source is needed.

This easy-to-use scanner fits very well with the Neoss Esthetic Healing Abutments and is what your practice needs. The best of all is that the NeoScan 1000 comes at the right price!

- ✓ Easy to use
- ✓ Fast scanning
- ✓ High precision
- ✓ Attractive price



“The NeoScan 1000 is a *superfast, lightweight, and easy-to-use* scanner. I had the pleasure of being part of early testing and have used the scanner for several digital impression indications at my clinic with excellent results. Digital dentistry is in need of more *cost-efficient* solutions so that clinicians can use it to its full potential. The NeoScan 1000 has the potential to do just this.”

– Dr. Marcus Dagnelid
Specialist in prosthodontics | Gothenburg, Sweden



Our offices

ARC SOLUTIONS

Arc Solutions
Box 13011
SE-250 13 Helsingborg
Sweden
T: +46 42 301 74 40
E: info@arcsolutions.se

AUSTRALIA & NEW ZEALAND

Neoss Australia Pty. Ltd
PO Box 404
New Farm
QLD 4005
Phone: +61 7 3216 0165
E: info.au@neoss.com

CHINA

Neoss China
7F, #303, Songhu Rd. (KIC Build #11)
Yangpu District
Shanghai, 200433
PR China
T: +86 21 80143135
E: china@neoss.com

EUROPEAN DISTRIBUTION

Neoss AG
The Circle 23
8058 Zurich-Airport
Switzerland
E: distributors@neoss.com

GERMANY & AUSTRIA

Neoss GmbH
Im Mediapark 5b
50670 Köln
T: +49 221 96980 10
F: +49 221 96980 199
E: info@neoss.de

ITALY

Neoss Italia S.r.l.
Viale Certosa 138
IT-20156 Milano
T: +39 02 92952 1
E: italia@neoss.info

MIDDLE EAST & AFRICA Neoss Limited (DMCC BRANCH)

Unit 1505 Fortune Tower
Cluster C – JLT
Dubai
P.O. Box 64093
United Arab Emirates
E: info.mea@neoss.com

SWEDEN, DENMARK & NORWAY

Neoss AB
Arvid Wallgrens backe 20
SE-413 46 Göteborg
T: +46 31 88 12 80
E: info@neoss.se

UNITED KINGDOM & REPUBLIC OF IRELAND

Neoss Ltd
Windsor House
Cornwall Road
Harrrogate
HG1 2PW
T: +44 1423 817 733
F: +44 1423 817 744
E: info@neoss.com

UNITED STATES

Neoss Inc.
21860 Burbank Blvd.,
Ste. 190 (North Entrance)
Woodland Hills, CA 91367
T: +1 866 626 3677
F: +1 818 432 2640
E: contact.usa@neoss.com

Our distributors

BAHRAIN

Cigalah Gulf Medical
Building 670, Block 309, Road 919
Salmuniya, Bahrain
T: +973 17232170
E: support@cgmed.com
maha@cgmed.com

CROATIA

Medical Intertrade d.o.o
10431 Sveta Nedelja
HRVATSKA
T: +385 1 33 74 020
F: +385 1 33 36 035
E: dental@medical-intertrade.hr

EGYPT

Dental Future
Building No 5
Mohamed Hassanein Heikal
St. Off Abbas El Akkad Street
1st Floor Flat 3
Cairo, Egypt
T: +20 109 6696449
E: mahmod@dentalfuturekw.com

IRAQ

Eblal Scientific Bureau
Al-Ameerat street
Almansour
Baghdad
Iraq
T: +964 780 411 0110
E: info@eblal.com

JAPAN

Neoss Japan Co.
Y101 Building 5F
No. 1-4 Kanda-Sudacho Chiyodaku
Tokyo 101-0041
T: +81 3 5289 3511
F: +81 3 5289 3512
E: info@neoss.co.jp

JORDAN

**Advanced Intelligence for Medical
Supplies & Equipment Ltd.**
Al-Swefeyyeh Paris str
Alhusainy Complex
Amman 11196
Jordan
T: +962 7 97389535
F: +962 6 5826797
E: f.khaddam@advanceime.com

KINGDOM OF SAUDI ARABIA

Tamer Group
Mohammad Ibn Abdul Wahab Street
Al Ammariya District
P.O. Box 180
Jeddah, 21411
Kingdom of Saudi Arabia
T: +966 12 643 5600
F: +966 12 642 4141
E: alaa.wanly@tamergruop.com

KUWAIT

Dental Future Company
Salmiya Salem el Mubark street
Suhaib complex office 10
Kuwait City
Kuwait
T: +965 51155997
E: Info@dentalfuturekw.com

LEBANON

Dentaltech Sarl
Brasilia Street
Hikmeh Building
First Basement
Baabca, Beirut – Lebanon
T: +961 5 950 707
F: +961 5 950 808

NETHERLANDS

Rotor, Medical Supplies B.V.
Rijksweg Zuid 51
NL-6071HV Swalmen
T: +31 46 481 12 02
F: +31 46 481 12 04
E: mail@rotor-medical.com

OMAN

Al Zahrawi Medical Supplies LLC.
Office No: 16 , Al Dar building opposite
Al Zubair automotive, street 61, Ghala
area.
Muscat, Oman
P.O. Box: 1627
T: +968 24587308
F: +968 24587487
E: info.oman@zahrawigroup.com

QATAR

Al Zahrawi Medical LLC
Building No. 36
Zone No. 39
Al Sadd Area
Street No.808
Doha, Qatar, P.O. Box 63065
T: +974 4431 872
F: +974 44365216
E: info.qatar@zahrawigroup.com

SPAIN

Formación, Gestión Y Rendimiento, S.L.
C. Nil Fabra 34-36.
Etslo 2 B
08012. Barcelona
Spain
T: +34 931 408 770
E: info@fgdental.com

SWITZERLAND

Flexident AG
Amlehstr. 22
CH-6010 Kriens
T: +41 41 3104020
E: info@flexident.ch

THAILAND

S.D. Tantavade Co., Ltd.
47 Samakkee 20, Samakkee Road
Thasai Muang
Nonthaburi 11000
Thailand
T: +66 2 952 4815 18
F: +66 2 952 4819
E: info@sdt1988.com

TURKEY

**Dentmark Sağlık Ürünleri
Anonim Şirketi**
Cevizli Mah. Cevizli D-100 Güney
Yanyol Cad.
No:2 Kuriş Kule Kat:3 D:52
Kartal-Istanbul PO: 34865
T: +90 216 417 4040
E: info@dentmark.com.tr

UNITED ARAB EMIRATES (ABU DHABI)

Elwan Technical Supplies
Airport Street Abu Dhabi UAE
T: +971 26393292
E: info@elwandental.com

UNITED ARAB EMIRATES (DUBAI & NORTHERN EMIRATES)

Al Zahrawi Medical Supplies
1st Floor, Ford Middle East &
Africa Building
Al Wasl Office Complex, Sheikh Zayed
Road, Downtown
Dubai, United Arab Emirates
P.O. Box: 5973
T: +971 4 262 2728
F: +971 4 262 5506
E: info@zahrawigroup.com

Contact your local representative for product availability details.
All products are not available in every market.