



Product Catalog

2025 | 2026

Intelligent Simplicity

Explore Intelligent Simplicity

Neoss offers intelligent solutions that are intuitively easy to use. Our products allow dental professionals to provide reliable and cost-effective treatments to their patients with predictable long-term results. With ingenuity and ambition, we strive to set new standards, leading the market with integrity. In developing smart treatment solutions and working closely with our customers, Neoss makes the complex less complicated. We call that Intelligent Simplicity.



Implants

- 4 Neoss ProActive® Edge Implants
- 6 Neoss ProActive® Tapered Implants
- Neoss ProActive® Wide Implants
- Neoss ProActive® Sinus Implants
- 7 Neoss ProActive® Straight Implants
- 8 Neoss ProActive® Straight NP Ø3.25 Implants
- Cover Screw
- Healing Abutments
- Esthetic Healing Abutments with ScanPeg

Instruments & More

- 9 Drills
- 11 Surgical Kits
- Surgical & Prosthetic Trays
- Surgical & Prosthetic Kits with instruments and drills
- 12 Neoss Guide Kit
- Neoss Guide Kit – SPARE PARTS
- Motors

Biomaterials

- 13 Instruments
- 14 Bone Allograft
- Mineralized Cortical/Cancellous Particulate
- Deminerlized bone fibers and mineralized bone particulate with Cortical/Cancellous Particulate **NEW**
- Bone Xenograft
- NeoGenix XP in jar
- NeoGenix XP in syringes
- 15 Resorbable Membranes
- Collagen Membranes
- Wound Dressings
- Non-Resorbable Membranes
- Ti-reinforced PTFE Membranes
- NeoGen® Cape
- 16 Non-Reinforced PTFE Membranes
- Fixation Kit
- Sutures

Prosthetics

- 17 Neoss4+ Treatment Solution
- Multi-Unit Abutments
- 18 Access Abutments
- 19 Multi-Unit/Access/Wide Screws **NEW**
- Provisional Abutments
- Esthetic Tissue Formers
- Provisional Titanium Abutments
- 20 Final Abutments
- Prepable Titanium Abutments
- Prepable Titanium Abutments, Alternative Emergence Profiles
- 21 Zirconia Abutments
- Final Restorations – Casting
- Gold NeoLink® (Direct Casting)
- CoCr Abutments (Indirect Casting)
- Burnout Abutments (Direct Casting)
- 22 Overdenture Abutments
- Equator Abutments
- Bar Abutment
- 23 Overdenture LOCATOR® Solution
- Impression Copings & Replicas
- Impression Copings
- 24 Replicas
- Screws **NEW**

Digital


- 25 Resonance Frequency Device (RFD)
- NeoTell™
- Intraoral Scanners
- NeoScan™ 1000
- NeoScan™ 2000
- 26 Digital Impressions
- Esthetic Healing Abutments with ScanPeg
- ScanPeg
- Scan Bodies **NEW**
- Scan Keys **NEW**
- Digital Analogs **NEW**
- 27 Digital Prosthetics
- NeoBase® Abutments SSC (Straight Screw Channel)
- NeoBase® Abutments ASC (Angulated Screw Channel)
- Neoss TiBase Abutments (compatible with Dentsply Sirona®)
- Ti NeoLink®

Support

- 29 Prosthetic Solution Standard and Narrow Platform, Overview
- Torque and Speed Recommendations, Overview
- 30 Symbols, Names and Abbreviations
- You can now buy your products online!
- Disclaimer of Liability


Neoss ProActive® implants – *Proven predictability* in implant treatment since 2002

All implants produced of highest quality, Commercially Pure Titanium Grade IV.



Neoss Implant

The Implant includes Implant and Cover Screw.




Neoss ProActive® Edge Implants

Art. No.	Description	A	B
21417	Implant – Neoss ProActive® Edge Ø3.5×9 mm		
21418	Implant – Neoss ProActive® Edge Ø3.5×11 mm		
21419	Implant – Neoss ProActive® Edge Ø3.5×13 mm	4.0 mm	1.6 mm
21420	Implant – Neoss ProActive® Edge Ø3.5×15 mm		
21421	Implant – Neoss ProActive® Edge Ø3.5×17 mm		
21423	Implant – Neoss ProActive® Edge Ø4.0×9 mm		
21424	Implant – Neoss ProActive® Edge Ø4.0×11 mm		
21425	Implant – Neoss ProActive® Edge Ø4.0×13 mm	4.3 mm	1.9 mm
21426	Implant – Neoss ProActive® Edge Ø4.0×15 mm		
21427	Implant – Neoss ProActive® Edge Ø4.0×17 mm		
21429	Implant – Neoss ProActive® Edge Ø4.5×9 mm		
21430	Implant – Neoss ProActive® Edge Ø4.5×11 mm		
21431	Implant – Neoss ProActive® Edge Ø4.5×13 mm	4.9 mm	2.2 mm
21432	Implant – Neoss ProActive® Edge Ø4.5×15 mm		
21433	Implant – Neoss ProActive® Edge Ø4.5×17 mm		
21435	Implant – Neoss ProActive® Edge Ø5.0×9 mm		
21436	Implant – Neoss ProActive® Edge Ø5.0×11 mm		
21437	Implant – Neoss ProActive® Edge Ø5.0×13 mm	5.4 mm	2.4 mm
21438	Implant – Neoss ProActive® Edge Ø5.0×15 mm		
21439	Implant – Neoss ProActive® Edge Ø5.0×17 mm		



SP = Standard Platform
Ø3.5 – 6.5 Implants





Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments in PEEK in two different sizes and Healing Abutment Screw.

Neoss ProActive® Edge Implants


Art. No.	Description	A	B
21301	Implant Kit – Neoss ProActive® Edge Ø3.5×9mm		
21302	Implant Kit – Neoss ProActive® Edge Ø3.5×11mm	4.0mm	1.6mm
21303	Implant Kit – Neoss ProActive® Edge Ø3.5×13mm		
21304	Implant Kit – Neoss ProActive® Edge Ø3.5×15mm		
21305	Implant Kit – Neoss ProActive® Edge Ø4.0×9mm		
21306	Implant Kit – Neoss ProActive® Edge Ø4.0×11mm	4.3mm	1.9mm
21307	Implant Kit – Neoss ProActive® Edge Ø4.0×13mm		
21308	Implant Kit – Neoss ProActive® Edge Ø4.0×15mm		
21309	Implant Kit – Neoss ProActive® Edge Ø4.5×9mm		
21310	Implant Kit – Neoss ProActive® Edge Ø4.5×11mm	4.9mm	2.2mm
21311	Implant Kit – Neoss ProActive® Edge Ø4.5×13mm		
21312	Implant Kit – Neoss ProActive® Edge Ø4.5×15mm		
21313	Implant Kit – Neoss ProActive® Edge Ø5.0×9mm		
21314	Implant Kit – Neoss ProActive® Edge Ø5.0×11mm	5.4mm	2.4mm
21315	Implant Kit – Neoss ProActive® Edge Ø5.0×13mm		
21316	Implant Kit – Neoss ProActive® Edge Ø5.0×15mm		





Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments (Ø5.0 – 6.5 implant kits include a wider version) in PEEK in two different sizes and Healing Abutment Screw.



Neoss ProActive® Tapered Implants

Art. No.	Description	A	B
21221	Implant Kit – Neoss ProActive® Tapered Ø3.5×9 mm		
21222	Implant Kit – Neoss ProActive® Tapered Ø3.5×11 mm	4.0 mm	2.1 mm
21223	Implant Kit – Neoss ProActive® Tapered Ø3.5×13 mm		
21224	Implant Kit – Neoss ProActive® Tapered Ø3.5×15 mm		
21227	Implant Kit – Neoss ProActive® Tapered Ø4.0×9 mm		
21228	Implant Kit – Neoss ProActive® Tapered Ø4.0×11 mm	4.3 mm	2.3 mm
21229	Implant Kit – Neoss ProActive® Tapered Ø4.0×13 mm		
21230	Implant Kit – Neoss ProActive® Tapered Ø4.0×15 mm		
21233	Implant Kit – Neoss ProActive® Tapered Ø4.5×9 mm		
21234	Implant Kit – Neoss ProActive® Tapered Ø4.5×11 mm	4.9 mm	2.8 mm
21235	Implant Kit – Neoss ProActive® Tapered Ø4.5×13 mm		
21236	Implant Kit – Neoss ProActive® Tapered Ø4.5×15 mm		
21239	Implant Kit – Neoss ProActive® Tapered Ø5.0×9 mm		
21240	Implant Kit – Neoss ProActive® Tapered Ø5.0×11 mm	5.4 mm	3.3 mm
21241	Implant Kit – Neoss ProActive® Tapered Ø5.0×13 mm		
21242	Implant Kit – Neoss ProActive® Tapered Ø5.0×15 mm		
21245	Implant Kit – Neoss ProActive® Tapered Ø5.5×9 mm		
21246	Implant Kit – Neoss ProActive® Tapered Ø5.5×11 mm	5.9 mm	3.8 mm
21247	Implant Kit – Neoss ProActive® Tapered Ø5.5×13 mm		

Neoss ProActive® Wide Implants


Art. No.	Description	A	B
21252	Implant Kit – Neoss ProActive® Ø6.0×7 mm		4.5 mm
21250	Implant Kit – Neoss ProActive® Ø6.0×9 mm	6.0 mm	4.3 mm
21251	Implant Kit – Neoss ProActive® Ø6.0×11 mm		4.3 mm

Neoss ProActive® Sinus Implants

Art. No.	Description	A
21211	Implant Kit – Neoss ProActive® Ø6.5×7 mm	
21212	Implant Kit – Neoss ProActive® Ø6.5×9 mm	6.7 mm
21213	Implant Kit – Neoss ProActive® Ø6.5×11 mm	



SP = Standard Platform
Ø3.5 – 6.5 Implants



Neoss Implant Kit

The Implant Kit includes Implant, Cover Screw, Healing Abutments (Ø5.0 – 5.5 implant kits include a wider version) in PEEK in two different sizes and Healing Abutment Screw.

Neoss ProActive® Straight Implants

Art. No.	Description	A	B
21181	Implant Kit – Neoss ProActive® Straight Ø3.5×7 mm	4.0 mm	2.4 mm
21182	Implant Kit – Neoss ProActive® Straight Ø3.5×9 mm		
21183	Implant Kit – Neoss ProActive® Straight Ø3.5×11 mm		
21184	Implant Kit – Neoss ProActive® Straight Ø3.5×13 mm		
21185	Implant Kit – Neoss ProActive® Straight Ø3.5×15 mm		
21186	Implant Kit – Neoss ProActive® Straight Ø3.5×17 mm		
21187	Implant Kit – Neoss ProActive® Straight Ø4.0×7 mm	4.0 mm	2.8 mm
21188	Implant Kit – Neoss ProActive® Straight Ø4.0×9 mm		
21189	Implant Kit – Neoss ProActive® Straight Ø4.0×11 mm		
21190	Implant Kit – Neoss ProActive® Straight Ø4.0×13 mm		
21191	Implant Kit – Neoss ProActive® Straight Ø4.0×15 mm		
21192	Implant Kit – Neoss ProActive® Straight Ø4.0×17 mm		
21193	Implant Kit – Neoss ProActive® Straight Ø4.5×7 mm	4.5 mm	3.2 mm
21194	Implant Kit – Neoss ProActive® Straight Ø4.5×9 mm		
21195	Implant Kit – Neoss ProActive® Straight Ø4.5×11 mm		
21196	Implant Kit – Neoss ProActive® Straight Ø4.5×13 mm		
21197	Implant Kit – Neoss ProActive® Straight Ø4.5×15 mm		
21198	Implant Kit – Neoss ProActive® Straight Ø4.5×17 mm		
21199	Implant Kit – Neoss ProActive® Straight Ø5.0×7 mm	5.0 mm	3.7 mm
21200	Implant Kit – Neoss ProActive® Straight Ø5.0×9 mm		
21201	Implant Kit – Neoss ProActive® Straight Ø5.0×11 mm		
21202	Implant Kit – Neoss ProActive® Straight Ø5.0×13 mm		
21203	Implant Kit – Neoss ProActive® Straight Ø5.0×15 mm		
21205	Implant Kit – Neoss ProActive® Straight Ø5.5×7 mm	5.5 mm	4.2 mm
21206	Implant Kit – Neoss ProActive® Straight Ø5.5×9 mm		
21207	Implant Kit – Neoss ProActive® Straight Ø5.5×11 mm		
21208	Implant Kit – Neoss ProActive® Straight Ø5.5×13 mm		




Implants & Healing Abutments



Neoss Implant Kit

The Neoss ProActive® NP Ø3.25 implant kit includes Implant, Cover Screw, Healing Abutment and Healing Abutment Screw.




Neoss ProActive® Straight NP Ø3.25 Implants

Art. No.	Description	A	B
21176	Implant Kit – Neoss ProActive® Straight, NP Ø3.25×9 mm		
21177	Implant Kit – Neoss ProActive® Straight, NP Ø3.25×11 mm	3.5mm	2.2mm
21178	Implant Kit – Neoss ProActive® Straight, NP Ø3.25×13 mm		
21179	Implant Kit – Neoss ProActive® Straight, NP Ø3.25×15 mm		

Cover Screws

Art. No.	Description
31408	Cover Screw SP
31409	Cover Screw NP



Healing Abutments

Healing Abutments


Art. No.	Description	A	B
31159	Healing Abutment Ti 2mm	2mm	
31160	Healing Abutment Ti 4mm	4mm	
31161	Healing Abutment Ti 6mm	6mm	4mm
31162	Healing Abutment Ti 8mm	8mm	
31163	Healing Abutment Ti 10mm	10mm	
31410	Healing Abutment PEEK 2.7 mm with Screw	2.7mm	5mm
31148	Healing Abutment PEEK 5mm with Screw	5mm	5mm
31320	Healing Abutment PEEK 5mm with Screw, NP		4mm
31321	Wide Healing Abutment PEEK 3.3mm with Screw	3.3mm	6.5mm
31335	Wide Healing Abutment PEEK 5mm with Screw	5mm	6.5mm

Esthetic Healing Abutments with ScanPeg

All Esthetic Healing Abutments include Titanium Screw and a ScanPeg

Art. No.	Description	A	B	C
31360	Esthetic Healing Abutment with ScanPeg – Wide Incisor		5.7mm	6.1mm
31361	Esthetic Healing Abutment with ScanPeg – Narrow Incisor		5.2mm	5mm
31362	Esthetic Healing Abutment with ScanPeg – Canine	5mm	6.3mm	5.4mm
31363	Esthetic Healing Abutment with ScanPeg – Premolar		6.8mm	5.8mm
31364	Esthetic Healing Abutment with ScanPeg – Molar		7.8mm	7.3mm

The instruments *you need*



Drills


Art. No. Description

41153	Round Bur Ø1.8 mm
41152	Lance Drill Ø 1.5×3/5/7 mm
41170	Pilot Drill Ø2.2/3.0 mm
41142	Drill Stop 7 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41143	Drill Stop 9 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41144	Drill Stop 11 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41145	Drill Stop 13 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41146	Drill Stop 15 mm Ø2.2, 3.0, 3.4, 3.9 and 4.4 mm
41205	Twist Drill, Straight Ø6.0×7 – 11 (13) mm
41215	Twist Drill, Straight Ø2.2×7 – 17 mm
41216	Twist Drill, Straight Ø2.85×7 – 17 mm
41217	Twist Drill, Straight Ø3.0×7 – 17 mm
41218	Twist Drill, Straight Ø3.2×7 – 17 mm
41219	Twist Drill, Straight Ø3.4×7 – 17 mm
41220	Twist Drill, Straight Ø3.6×7 – 17 mm
41221	Twist Drill, Straight Ø3.9×7 – 17 mm
41222	Twist Drill, Straight Ø4.1×7 – 17 mm
41223	Twist Drill, Straight Ø4.4×7 – 17 mm
41224	Twist Drill, Straight Ø4.6×7 – 17 mm
41225	Twist Drill, Straight Ø4.9×7 – 17 mm
41226	Twist Drill, Straight Ø5.1×7 – 17 mm

Art. No. Description


41233	Drill Kit, Straight Implants Ø3.5 – 5.0 <i>Includes Twist Drills Ø2.2, 3.0, 3.4, 3.9 and 4.4</i>
41184	Drill Kit, Straight – Short 7 – 13 mm <i>Includes Twist Drills Ø2.2, 3.0, 3.4, 3.9, 4.4 and 4.9</i>
41227	Twist Drill, Tapered Ø3.0×7 – 17 mm
41228	Twist Drill, Tapered Ø3.4×7 – 17 mm
41229	Twist Drill, Tapered Ø3.9×7 – 17 mm
41230	Twist Drill, Tapered Ø4.4×7 – 17 mm
41231	Twist Drill, Tapered Ø4.9×7 – 17 mm
41232	Twist Drill, Tapered Ø5.5×7 – 17 mm
41235	Drill Kit, Tapered Implants Ø3.5 – 5.0 <i>Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T and 4.4T</i>
41191	Drill Kit Tapered – Short 7 – 13 mm <i>Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T, 4.4T, and 4.9T</i>
41236	Drill Kit, Edge Implants Ø3.5 – 5.0 <i>Includes Twist Drills Ø2.2S, 3.0T, 3.4T, 3.9T and 4.4T</i>

41139	Countersink, Straight Ø3.25	2.7mm
41110	Countersink, Straight Ø3.5	2.9mm
41111	Countersink, Straight Ø4.0	3.3mm
41112	Countersink, Straight Ø4.5	3.5mm
41133	Countersink, Straight Ø5.0	4.3mm
41127	Countersink, Straight Ø5.5	4.7mm
41178	Countersink, Tapered Ø3.5	2.0mm
41179	Countersink, Tapered Ø4.0	2.8mm
41180	Countersink, Tapered Ø4.5	3.3mm
41181	Countersink, Tapered Ø5.0	3.7mm
41182	Countersink, Tapered Ø5.5	4.2mm
41210	Countersink, Edge Ø3.5	1.8mm
41211	Countersink, Edge Ø4.0	2.6mm
41212	Countersink, Edge Ø4.5	3.0mm
41213	Countersink, Edge Ø5.0	3.5mm
41202	Countersink Ø6.0	4.8mm
41201	Countersink Ø6.5	5.2mm
41113	Screw Tap Ø3.5	
41114	Screw Tap Ø4.0	
41115	Screw Tap Ø4.5	
41134	Screw Tap Ø5.0	
41126	Screw Tap Ø5.5	
41140	Screw Tap Ø3.25	
41237	Bone Profiler Ø5.1mm, SP	
41238	Bone Profiler Ø6.6mm, SP	
41239	Bone Profiler Ø4.1 mm, NP	
41240	Bone Remover Cover Screw, SP	
41241	Bone Remover Cover Screw, NP	




Neoss Prosthetic Tray
– NeoNeoss ProActive*
Implant System Box

ProActive* Edge



ProActive* Tapered



ProActive* Straight



Neoss Prosthetic Instrument Kit – Neo



Drill and Surgical Instrument Kit, Edge



Drill and Surgical Instrument Kit, Tapered



Drill and Surgical Instrument Kit, Straight

Surgical Kits

Surgical & Prosthetic Trays

Art. No. Description

51180	Neoss Prosthetic Tray – Neo
51186	Neoss ProActive* Implant System Box
51189	Neoss Surgical and Prosthetic Tray – ProActive* Edge
51188	Neoss Surgical and Prosthetic Tray – ProActive* Tapered
51187	Neoss Surgical and Prosthetic Tray – ProActive* Straight

Surgical & Prosthetic Kits with instruments and drills

Art. No. Description

51181	Neoss Prosthetic Instrument Kit – Neo <i>Includes 51180, 51121, 51139, 51140, 51141 and 51126</i>
41151	Surgical Instrument Kit <i>Includes art no 41120, 51118, 51119, 51144, 51139, 51140, 51125, 51126</i>
41253	Drill and Surgical Instrument Kit, Edge Implants Ø3.5 – 5.0 <i>Multiple use twist drills. Includes art no 41236, 41210, 41211, 41212, 41213, 41237, 41238, 41240, 41151, 41153, 41152, 41170, 51121, 51137, 51141, 51186, 51189 and 51190</i>
41251	Drill and Surgical Instrument Kit, Tapered Implants Ø3.5 – 5.0 <i>Multiple use twist drills. Includes art no 41235, 41218, 41220, 41222, 41224, 41237, 41238, 41240, 41178, 41179, 41180, 41181, 41151, 41153, 41152, 41170, 51121, 51137, 51141, 51186, 51188 and 51190</i>
41247	Drill and Surgical Instrument Kit, Straight Implants Ø3.5 – 5.0 <i>Multiple use twist drills. Includes art no 41233, 41218, 41220, 41222, 41224, 41110, 41111, 41112, 41133, 41151, 41170, 41153, 41152, 41237, 41238, 41240, 51121, 51137, 51141, 51186, 51187 and 51190</i>
41249	Drill and Surgical Instrument Kit, Straight Implants NP Ø3.25 <i>Multiple use twist drills. Includes art no 41215, 41216, 41217, 41139, 51142, 51145, 41241, 41239, 41170, 41153, 41152, 51121, 51137, 41120, 51144, 51125, 51139, 51140, 51141, 51126, 51190, 51186 and 51187</i>
41255	Drill and Surgical Instrument Add-On-Kit, Straight Implant Ø3.25 <i>Multiple use twist drills. Includes art no 41216, 41239, 41241, 41139, 51142 and 51145</i>



Drill and Surgical Instrument Kits also includes the prosthetic instruments and a Depth Gauge Probe. All in one box for your convenience.

Neoss® Guide Kit

Art. No. Description

51154	Neoss® Guide Kit Ø3.5 – Ø6.0 (Neoss Implants Ø3.5 – Ø6.0mm, 7 – 15mm lengths.) Includes 5 × Drill-Hubs, 2 × Countersink Conduits, 4 × Implant Mount Ø5.0 S9, 2 × Implant Mount Ø6.0 S9, 7 × Guide Keys, 3 × Guide Key Adapter, 1 × Holder Guide Key Adapter. Spare parts available, contact local representative for details
51173	Neoss® Guide – Implant Mount SP, Ø5.0 S9 (2 pcs)
51174	Neoss® Guide – Implant Mount SP, Ø6.0 S9 (2 pcs)
51156	Neoss® Guide – Sleeve Ø5.0 (10 pcs)
51157	Neoss® Guide – Sleeve Ø6.0 (5 pcs)
51155	Neoss® Guide – Implant Inserters, SP For partially guided only, see IFU




Neoss® Guide Kit



Implant Mount



Sleeve




Implant Inserters

Neoss® Guide Kit – SPARE PARTS

Art. No. Description

51160	Neoss® Guide – Drill-Hub I
51161	Neoss® Guide – Drill-Hub II
51162	Neoss® Guide – Drill-Hub III
51163	Neoss® Guide – Drill-Hub IV
51164	Neoss® Guide – Drill-Hub V
51165	Neoss® Guide – Guide Key Ø2.2/2.85
51166	Neoss® Guide – Guide Key Ø3.0/3.2
51167	Neoss® Guide – Guide Key Ø3.4/3.6
51168	Neoss® Guide – Guide Key Ø3.9/4.1
51169	Neoss® Guide – Guide Key Ø4.4/4.6
51170	Neoss® Guide – Guide Key Ø4.9/5.1
51178	Neoss® Guide – Guide Key Ø5.5
51171	Neoss® Guide – Countersink Conduit Ø5.0
51172	Neoss® Guide – Countersink Conduit Ø6.0
51176	Neoss® Guide – Guide Key Adapter (2 pcs)
51177	Neoss® Guide – Holder Guide Key Adapter
51158	Drill Sleeve Ø2.2 (10 pcs)



Drill-Hubs



Guide Key



Countersink Conduit



Drill Sleeve




Guide Key Adapter



Holder Guide Key Adapter

Motors

Neoss offers a diverse portfolio of high-quality motors designed to meet your specific needs. For more information, please don't hesitate to contact us. Our product range includes renowned manufacturers such as W&H and Saeshin, among others.



Instruments

Art. No. Description

51137	Implant Inserter, SP 17 mm
51118	Implant Inserter, SP 22 mm
51119	Implant Inserter, SP 32 mm
51145	Implant Inserter, NP 24 mm
51142	Implant Inserter, NP 32 mm
51139	Neo Screwdriver, Machine 22 mm
51140	Neo Screwdriver, Machine 32 mm
51141	Neo Screwdriver, Manual 15 mm
51183	iGO Screwdriver, Machine 22 mm
51184	iGO Screwdriver, Machine 32 mm
51125	Direction/Depth Gauge (4 pcs)
41120	Drill Extender
51190	Depth Gauge Probe
51121	Ratchet
51126	Manual Handle <i>To be used with prosthetic drivers only</i>
51191	Insertion Handle – Large
51144	Wrench Adapter <i>To be used with implant inserters only</i>
51149	Angulation Gauge, SP
90138	Adjustment Handle

NEO

iGO

Solutions for your guided bone *regeneration* procedures

Bone Allograft

Mineralized Cortical/Cancellous Particulate

Art. No.	Description
64063	NeoGenix™ Allograft Mineralized Cortical Cancellous Particulate 250-1000µm 0.25cc-U
64064	NeoGenix™ Allograft Mineralized Cortical Cancellous Particulate 250-1000µm 0.5cc-U
64065	NeoGenix™ Allograft Mineralized Cortical Cancellous Particulate 250-1000µm 1.0cc-U
64066	NeoGenix™ Allograft Mineralized Cortical Cancellous Particulate 250-1000µm 2.0cc-U
64067	NeoGenix™ Allograft Mineralized Cortical Cancellous Particulate 250-1000µm 2.5cc-U
64068	NeoGenix™ Allograft Mineralized Cortical Cancellous Particulate 250-1000µm 5.0cc-U



NeoGenix™ Allograft Mineralized

Demineralized bone fibers and mineralized bone particulate with Cortical/Cancellous Particulate

NEW

Art. No.	Description	Height	Diameter
64075	NeoGenix™ Extraction Allograft	11mm	7.5mm



NeoGenix™ Extraction Allograft

Bone Xenograft

NeoGenix XP™ in jar

Art. No.	Description	Weight	Volume
64091	NeoGenix XP™ 0.5 cc (0.25 – 1 mm particle size) Jar	0.15 g	0.5 cc
64092	NeoGenix XP™ 1.0 cc (0.25 – 1 mm particle size) Jar	0.25 g	1.0 cc
64093	NeoGenix XP™ 2.0 cc (0.25 – 1 mm particle size) Jar	0.50 g	2.0 cc
64094	NeoGenix XP™ 4.0 cc (0.25 – 1 mm particle size) Jar	1.00 g	4.0 cc
64095	NeoGenix XP™ 1.0 cc (1 – 2 mm particle size) Jar	2.00 g	1.0 cc
64096	NeoGenix XP™ 2.0 cc (1 – 2 mm particle size) Jar	0.50 g	2.0 cc



NeoGenix XP in jar

NeoGenix XP™ in syringes

Art. No.	Description	Weight	Volume
64097	NeoGenix XP™ 0.25 cc (0.25 – 1 mm particle size) Curved Syringe	0.10 g	0.25 cc
64098	NeoGenix XP™ 0.5 cc (0.25 – 1 mm particle size) Curved Syringe	0.21 g	0.5 cc



NeoGenix XP in syringes

Resorbable Membranes

Collagen Membranes

Art. No. Description

NeoGen® Collagen Flex resorbs in 3–4 months

		A	B
64050	NeoGen® Collagen Membrane Flex S	20 mm	15 mm
64051	NeoGen® Collagen Membrane Flex M	30 mm	20 mm
64052	NeoGen® Collagen Membrane Flex L	40 mm	30 mm

NeoGen® Collagen Firm resorbs in 6–9 months

		A	B
64055	NeoGen® Collagen Membrane Firm S	20 mm	15 mm
64056	NeoGen® Collagen Membrane Firm M	30 mm	20 mm
64057	NeoGen® Collagen Membrane Firm L	40 mm	30 mm



Collagen Flex

Collagen Firm






Wound Dressings

Wound Dressings

Art. No. Description

Wound dressings essentially resorb within 30 days of placement

		A	B
64060	NeoGen® Collagen Wound Dressing Tape (1 mm thick, 10 pcs)	75 mm	25 mm
64061	NeoGen® Collagen Wound Dressing Foam (4 mm thick, 10 pcs)	40 mm	20 mm
64062	NeoGen® Collagen Wound Dressing Plug (10 pcs)	20 mm	10 mm



Ti-Reinforced Membranes

Non-Resorbable Membranes

Ti-Reinforced PTFE Membranes

Art. No. Description

		A	B	C	D
64010	NeoGen® PTFE Membrane, Ti-Reinforced – S I (Small Interproximal)	29 mm	14 mm	19 mm	4 mm
64011	NeoGen® PTFE Membrane, Ti-Reinforced – M I (Medium Interproximal)	30 mm	19 mm	19 mm	9 mm
64012	NeoGen® PTFE Membrane, Ti-Reinforced – L I (Large Interproximal)	36 mm	21 mm	26 mm	11 mm
64013	NeoGen® PTFE Membrane, Ti-Reinforced – M (Medium)	32 mm	22 mm	23 mm	12 mm
64014	NeoGen® PTFE Membrane, Ti-Reinforced – L (Large)	34 mm	25 mm	24 mm	15 mm
64015	NeoGen® PTFE Membrane, Ti-Reinforced – XL (Extra Large)	46 mm	32 mm	24 mm	15 mm

NeoGen® Cape

All NeoGen® Cape Spacer include a NeoGen® Cape Screw

Art. No. Description

		A	B	C	D
64040	NeoGen® Cape PTFE Membrane – M	28 mm	30 mm	16 mm	9 mm
64041	NeoGen® Cape PTFE Membrane – L	28 mm	30 mm	22 mm	9 mm
64080	NeoGen® Cape Spacer 0.7 mm, Neoss SP				
64081	NeoGen® Cape Spacer 1.5 mm, Neoss SP				
64082	NeoGen® Cape Spacer 0.7 mm, Neoss NP				
64083	NeoGen® Cape Spacer 1.5 mm, Neoss NP				



NeoGen® Cape

NeoGen® Cape Spacer

NeoGen® Cape Screw


Non-Reinforced PTFE Membranes

Art. No.	Description	A	B
64019	NeoGen® PTFE Membrane, Non-Reinforced	18mm	12mm
64020	NeoGen® PTFE Membrane, Non-Reinforced	29mm	14mm
64021	NeoGen® PTFE Membrane, Non-Reinforced	34mm	25mm




Fixation Kit

Art. No.	Description
64031	Tack Positioning Instrument
64033	Tack Mallet
64034	Tack and Screw Cassette
64036	Membrane Tack 3mm (10 pcs)
64037	Membrane Screw 3mm (10 pcs)




Sutures

Art. No.	Description
32.Z4401.00	Neoss PTFE Suture 60 cm 4/0 – Reverse Cutting 16 mm 3/8 (12 pcs)
32.Z4408.00	Neoss PTFE Suture 60 cm 6/0 – Reverse Cutting 12 mm 3/8 (12 pcs)



Experience *easy-to-use full arch flexibility* with predictable results using Neoss4+ treatment solution



Neoss4+ treatment solution


Multi-Unit Abutments

*Clinical screw fixed within the Multi-Unit Abutment
Healing Cap and Replica needs to be ordered separately*

Art. No.	Description	A	B	C
31411	Multi-Unit Abutment 1.0 mm	1 mm	4.6 mm	
31412	Multi-Unit Abutment 2.0 mm	2 mm		
31413	Multi-Unit Abutment 3.0 mm	3 mm		
31414	Multi-Unit Abutment 4.0 mm	4 mm		
31425	Multi-Unit Abutment 10° 2.0 mm	2 mm	4.6 mm	
31426	Multi-Unit Abutment 10° 3.0 mm	3 mm		
31427	Multi-Unit Abutment 10° 4.0 mm	4 mm		
31429	Multi-Unit Abutment 17° 2.5 mm	2.5 mm	4.6 mm	
31430	Multi-Unit Abutment 17° 3.5 mm	3.5 mm		
31431	Multi-Unit Abutment 17° 4.5 mm	4.5 mm		
31433	Multi-Unit Abutment 30° 3.0 mm	3 mm	4.6 mm	
31434	Multi-Unit Abutment 30° 4.0 mm	4 mm		
31435	Multi-Unit Abutment 30° 5.0 mm	5 mm		
31440	Multi-Unit Healing Cap	5 mm	5 mm	
31442	Multi-Unit Impression Coping Open tray 6 mm			
31443	Multi-Unit Impression Coping Closed tray 10 mm			
31444	Multi-Unit Replica (5 pcs)			
31445	Multi-Unit Provisional Abutment Ti <i>Includes Provisional/Laboratory Screw</i>	15.7 mm	0.7 mm	5 mm
31446	Multi-Unit NeoBase® SSC G0.5 mm – H4.5 mm <i>Includes Provisional/Laboratory Screw and final screw</i>	5 mm	0.5 mm	5.1 mm
31447	Multi-Unit NeoBase® ASC G0.5 mm – H4.5 mm iGO <i>Includes Provisional/Laboratory Screw and the iGO screw (requires the iGO screwdriver, 51183 or 51184)</i>	5 mm	0.5 mm	5.1 mm
31458	Scan Body Ti Neoss® Multi-Unit – 1 pc			
31462	Scan Body Ti Neoss® Multi-Unit – 5 pcs			
31473	Analog for Printed Model Neoss® Multi-Unit			

See page 26 for ScanKey.








Access Abutments

Clinical screw fixed within the Access Abutment

Art. No.	Description	A	B
All Access Abutments include Abutment and a Healing Abutment			
31238	Access Abutment 1.5mm	1.5mm	4mm
31215	Access Abutment 3.0mm	3mm	
31216	Access Abutment 4.0mm	4mm	
31218	Access Abutment 6.0mm	6mm	
31219	Access Abutment 10° 2.6mm	2.6mm	4mm
31243	Access Abutment 10° 4.6mm	4.6mm	
31221	Access Abutment 20° 2.6mm	2.6mm	
31244	Access Abutment 20° 4.6mm	4.6mm	
31223	Access Abutment 30° 2.9mm	2.9mm	
31245	Access Abutment 30° 4.9mm	4.9mm	
31225	Access Impression Coping 8mm and Replica		
31366	Access Impression Coping 13mm and Replica <i>31225 and 31366 include Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>		
31453	Scan Body Ti Neoss® Access – 1 pc		
31461	Scan Body Ti Neoss® Access – 5 pcs		
31472	Digital Analog Neoss® Access – 1 pc		
90207	Access Equator Abutment Kit 0.5mm		
90208	Access Equator Abutment Kit 1.5mm <i>90207 and 90208 include Equator Abutment, Housing with processing Cap, Protector Disk, Standard, Strong, Soft and Extra Soft Retention caps</i>		



Art. No.	Description	A	B	C
31226	Access Ti NeoLink® Multi	4.3mm	0.3mm	4.5mm
31227	Access Gold NeoLink® Multi <i>31226 and 31227 include NeoLink®, Prosthetic Screw, Provisional/Laboratory Screw and a Plastic Coping</i>	4.7mm	0.7mm	4.5mm
31301	Access Provisional Ti Abutment Multi <i>Includes Provisional/Laboratory Screw and Straight Coping</i>	15.7mm	0.7mm	4.5mm
31266	Access Burnout Abutment Multi – 2pcs <i>Includes Straight Coping</i>			




Multi-Unit/Access/Wide Screws **NEW**

Overview abutment, prosthetic and laboratory screws

	Laboratory and Provisional prosthetics				Definitive prosthetics			
	SP	NP	Multi-Unit / Access	Multi-Unit / Access for Zr/Temp bridges	SP	NP	Multi-Unit / Access	Multi-Unit / Access for Zr/Temp bridges
Straight screw channel	 31149 Neo Provisional/ Laboratory Screw, SP	 31239 Neo Provisional/ Laboratory Screw, NP	 31257 Multi-Unit/ Access Neo Provisional/ Laboratory Screw	 31480 Multi-Unit/ Access Neo Provisional/ Laboratory Wide Screw	 31258 Neo Abutment Screw, SP	 31240 Neo Abutment Screw Ti, NP	 31229 Multi-Unit/ Access Neo Prosthetic Screw	 31482 Multi-Unit/ Access Neo Prosthetic Wide Screw
Angulated screw channel	 31350 iGO Provisional/ Laboratory Screw, SP	 31352 iGO Provisional/ Laboratory Screw, NP	 31351 Multi-Unit/ Access iGO Provisional/ Laboratory Screw	 31484 Multi-Unit/ Access iGO Provisional/ Laboratory Wide Screw	 31347 iGO Abutment Screw, SP	 31349 iGO Abutment Screw, NP	 31348 Multi-Unit/ Access iGO Prosthetic Screw	 31486 Multi-Unit/ Access iGO Prosthetic Wide Screw

For our full range of screws, please see page 24.



Provisional Abutments

Esthetic Tissue Formers

All Esthetic Tissue Formers include a metal base that allows for radiographic verification, and include a provisional screw.

Art. No.	Description	A	B	C	D
31340	Esthetic Tissue Former – Wide Incisor			5.7mm	6.1mm
31341	Esthetic Tissue Former – Narrow Incisor			5.2mm	5mm
31342	Esthetic Tissue Former – Canine	10mm	4mm	6.3mm	5.3mm
31343	Esthetic Tissue Former – Premolar			6.8mm	5.8mm
31344	Esthetic Tissue Former – Molar			7.8mm	7.3mm




Esthetic Tissue Formers

Provisional Titanium Abutments


All Provisional Abutments include Provisional Screw and Coping

Art. No.	Description	A	B	C
31267	Provisional Ti Abutment Mono			
31268	Provisional Ti Abutment Multi	15.7mm	0.7mm	4.5mm
31332	Provisional Ti Abutment Mono – Cement Retained			
31241	Provisional Ti Abutment, NP	10.7mm	0.7mm	3.8mm

Mono - Includes Provisional Titanium Abutment, Provisional Screw NP



Provisional Ti Abutment, NP



Final Abutments

Preparable Titanium Abutments

All Preparable Abutments include Neo Abutment Screw and Neo Laboratory Screw

Art. No.	Description	A	B	C	D	E
31280	Preparable Ti Abutment 1.5 mm 0° – Wide Incisor	8.5mm	2.5mm	1.5mm	5.8mm	6.5mm
31281	Preparable Ti Abutment 1.5 mm 15° – Wide Incisor					
31282	Preparable Ti Abutment 1mm 0° – Narrow Incisor	8mm	2mm	1mm	5mm	5.1mm
31283	Preparable Ti Abutment 1mm 15° – Narrow Incisor					
31284	Preparable Ti Abutment 3mm 0° – Narrow Incisor	10mm	4mm	3mm	5.1mm	5.1mm
31285	Preparable Ti Abutment 3mm 15° – Narrow Incisor					
31286	Preparable Ti Abutment 3mm 0° – Canine	10mm	4mm	3mm	6.7mm	5.4mm
31305	Preparable Ti Abutment 1.5 mm 0° – Premolar	7mm	2mm	1.5mm	6.9mm	5.7mm
31287	Preparable Ti Abutment 3mm 0° – Premolar	9mm	3mm	3mm	6.9mm	5.7mm
31306	Preparable Ti Abutment 1.5 mm 0° – Molar	7mm	2mm	1.5mm	7.8mm	7.1mm
31288	Preparable Ti Abutment 3mm 0° – Molar	9mm	3mm	3mm	7.8mm	7.1mm
31201	Preparable Ti Abutment 1mm 0°, NP	8mm	2mm	1mm	4mm	
31202	Preparable Ti Abutment 1mm 15°, NP					



Preparable Titanium Abutments, Alternative Emergence Profiles

Art. No.	Description	A	B	C	D	E
31317	Preparable Ti Abutment 2mm 0° – Wide Emergence <i>Includes Neo Abutment Screw, Neo Laboratory Screw, specific Cover Screw and specific Healing Abutment PEEK with Screw. Compatible with Implant diameters Ø5.0 – 6.5mm</i>	9mm	2mm	2mm	7.9mm	7.1mm
31318	Preparable Ti Abutment 1.5 mm 0° – Narrow Emergence	8.5mm	2mm	1.5mm	4.4mm	4.5mm
31319	Preparable Ti Abutment 1.5 mm 15° – Narrow Emergence <i>All Narrow Emergence Preparable Abutments include Neo Abutment Screw and Neo Laboratory Screw. Only use with Ø3.5 and Ø4.0 mm implants</i>					




Zirconia Abutments

All Zirconia Abutments include separate Zirconia Coping, Ti NeoLink® Mono, Neo Abutment Screw and Neo Laboratory Screw

Art. No.	Description	A	B	C	D	E
31290	Zirconia Abutment – Wide Incisor				5.8mm	6.4mm
31291	Zirconia Abutment – Narrow Incisor		2.5mm	1.5mm	5.4mm	5.1mm
31292	Zirconia Abutment – Canine	9mm			6.2mm	5.8mm
31293	Zirconia Abutment – Premolar		2mm	2mm	6.4mm	5.5mm
31294	Zirconia Abutment – Molar				7.2mm	6.7mm






Final Restorations – Casting

Gold NeoLink® (Direct Casting/Aka UCLA Abutment)

Gold NeoLinks include Neo Abutment Screw, Neo Laboratory Screw and Straight Copings

Art. No.	Description	A	B	C
31102	Gold NeoLink® Mono	4.0mm	0.4mm	4.5mm
31103	Gold NeoLink® Multi			3.8mm
31199	Gold NeoLink® Mono, NP <i>Includes Neo Abutment Screw Ti NP, Neo Laboratory Screw NP and a set of Preformed Plastic Copings</i>			



CoCr Abutment Mono



CoCr Abutment Multi



Burnout Abutment Mono



Burnout Abutment Multi



Access Burnout Abutment Multi



CoCr Abutments (Indirect Casting)

CoCr Abutments include Neo Abutment Screw, Neo Laboratory Screw and Straight Coping

Art. No.	Description	A	B	C
31269	CoCr Abutment Mono	15.7mm	0.7mm	4.5mm
31260	CoCr Abutment Multi			

Burnout Abutments (Direct Casting)

31333, 31157 and 31266 include Straight Coping

Art. No.	Description	A	B	C
31333	Burnout Abutment Mono (2pcs)	15.7mm	0.7mm	4.5mm
31157	Burnout Abutment Multi (2pcs)			
31266	Access Burnout Abutment Multi (2pcs)			
31255	Burnout NeoLink® Multi NP (2pcs) <i>Includes set of Preformed Plastic Copings</i>	4.3mm		3.8mm

Overdenture Abutments

Equator Abutments

All kits include Equator Abutment, Housing with processing Cap, Protector Disk, Standard, Strong, Soft and Extra Soft Retention caps

Art. No. Description

90201	Neoss Equator Abutment Kit 1mm
90202	Neoss Equator Abutment Kit 2mm
90242	Neoss Equator Abutment Kit 2mm, NP
90203	Neoss Equator Abutment Kit 3mm
90243	Neoss Equator Abutment Kit 3mm, NP
90204	Neoss Equator Abutment Kit 4mm
90244	Neoss Equator Abutment Kit 4mm, NP
90205	Neoss Equator Abutment Kit 5mm
90206	Neoss Equator Abutment Kit 6mm
90207	Access Equator Abutment Kit 0.5mm
90208	Access Equator Abutment Kit 1.5mm
90210	Neoss Equator Cap and Housing Kit
90215	Neoss Equator Cap and Housing Kit – Extended Divergence
90211	Neoss Equator Cap Standard Retention, White (2pcs)
90212	Neoss Equator Cap Strong Retention, Violet (2pcs)
90213	Neoss Equator Cap Soft Retention, Pink (2pcs)
90214	Neoss Equator Cap Extra Soft Retention, Yellow (2pcs)
90220	Neoss Equator Impression Coping & Replica (2 sets)
90222	Neoss Equator Protector Disk (10 pcs)
90232	Neoss Equator Cap Tool



Bar Abutment

Art. No. Description

31168	Bar Abutment Gold
<i>Includes Abutment, Neo Abutment Screw and Neo Laboratory Screw</i>	





LOCATOR® Core Tool



LOCATOR® Impression Coping & Replica



LOCATOR® Replacement Males



LOCATOR® Block Out Spacer



LOCATOR® and Access LOCATOR® Abutment




Equator Abutment kit



Impression Coping and Replica



Impression Coping Open Tray and Replica



Impression Coping and Replica, NP


Multi-Unit Impression Coping
Open tray

Multi-Unit Impression Coping
Closed tray


Access Impression Coping and Replica



Impression Coping Cap



Impression Coping Extension Tube

Overdenture LOCATOR® Solution

Art. No. Description

08260	LOCATOR® 35 Ncm Torque Wrench Inserter Driver 15 mm
08393	LOCATOR® core tool
08451	1.0 mm LOCATOR® Abutment for Neoss 4.1 mm D Implant
08452	2.0 mm LOCATOR® Abutment for Neoss 4.1 mm D Implant
08453	3.0 mm LOCATOR® Abutment for Neoss 4.1 mm D Implant
08454	4.0 mm LOCATOR® Abutment for Neoss 4.1 mm D Implant
08455	5.0 mm LOCATOR® Abutment for Neoss 4.1 mm D Implant
08456	6.0 mm LOCATOR® Abutment for Neoss 4.1 mm D Implant
08505	LOCATOR® impression coping (4 pcs)
08519-2	LOCATOR® male processing package assembly: (2 pcs) <i>Package includes, block out spacer, housing, clear, pink and blue replacement males</i>
08524	LOCATOR® replacement male Clear (4 pcs)
08527	LOCATOR® replacement male Light Retention, Pink (4 pcs)
08529	LOCATOR® replacement male Extra Light Retention, Blue (4 pcs)
08530	LOCATOR® 4 mm female analog (4 pcs)
08547	LOCATOR® replacement male, Green (4 pcs) <i>Extended Range</i>
08548	LOCATOR® replacement male, Red (4 pcs) <i>Extra Light Retention Extended Range</i>
08913	LOCATOR® Driver Machine, 23 mm
90107	Access LOCATOR® Abutment Kit 0.5 mm

Impression Copings & Replicas

Impression Copings

Art. No. Description


31170	Impression Coping 8mm and Replica
31145	Impression Coping 11mm and Replica
31183	Impression Coping 18 mm and Replica <i>31170, 31145 and 31183 include Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>
31302	Impression Coping Open Tray 8mm and Replica
31185	Impression Coping Open Tray 13 mm and Replica <i>31302 and 31185 include Coping, Screw, Replica and Extension Tube</i>
31197	Impression Coping 11mm and Replica, NP <i>Includes Coping, Screw, Replica, Extension Tube for Open Tray and Plastic Cap for Closed Tray Technique</i>
31442	Multi-Unit Impression Coping Open tray 6mm*
31443	Multi-Unit Impression Coping Closed tray 10mm*
31366	Access Impression Coping 13mm and Replica
31186	Impression Coping Cap Closed Tray (10pcs)
31262	Impression Coping Extension Tube (20 pcs)

*Analog not included, should be ordered separately

Replicas

Art. No. Description

31444	Multi-Unit Replica (5 pcs)*
31231	Access Replica (5pcs)
31146	Implant Replica (5pcs)
31198	Implant Replica, NP (5pcs)
31304	Implant Replica, NP (1pc)
31166	Protection Replica (1pc)



Screws


Overview abutment, prosthetic and laboratory screws

	Laboratory and Provisional prosthetics				Definitive prosthetics			
	SP	NP	Multi-Unit / Access	Multi-Unit / Access for Zr/Temp bridges	SP	NP	Multi-Unit / Access	Multi-Unit / Access for Zr/Temp bridges
Straight screw channel 	 31149 Neo Provisional/Laboratory Screw, SP	 31239 Neo Provisional/Laboratory Screw, NP	 31257 Multi-Unit/Access Neo Provisional/Laboratory Screw	 31480 Multi-Unit/Access Neo Provisional/Laboratory Wide Screw	 31258 Neo Abutment Screw, SP	 31240 Neo Abutment Screw Ti, NP	 31229 Multi-Unit/Access Neo Prosthetic Screw	 31482 Multi-Unit/Access Neo Prosthetic Wide Screw
Angulated screw channel 	 31350 iGO Provisional/Laboratory Screw, SP	 31352 iGO Provisional/Laboratory Screw, NP	 31351 Multi-Unit/Access iGO Provisional/Laboratory Screw	 31484 Multi-Unit/Access iGO Provisional/Laboratory Wide Screw	 31347 iGO Abutment Screw, SP	 31349 iGO Abutment Screw, NP	 31348 Multi-Unit/Access iGO Prosthetic Screw	 31486 Multi-Unit/Access iGO Prosthetic Wide Screw

Art. No. Description

31258	Neo Abutment Screw, SP - 1pc
31259	Neo Abutment Screw, SP - 5pcs
31240	Neo Abutment Screw Ti, NP - 1pc
31149	Neo Provisional/Laboratory Screw, SP - 5pcs
31239	Neo Provisional/Laboratory Screw, NP - 1pc
31303	Neo Laboratory Screw Long, SP - 5pcs
31347	iGO Abutment Screw, SP - 1pc
31349	iGO Abutment Screw, NP - 1pc
31348	Multi-Unit/Access iGO Prosthetic Screw - 1pc
31350	iGO Provisional/Laboratory Screw, SP - 1pc
31352	iGO Provisional/Laboratory Screw, NP - 1pc
31351	Multi-Unit/Access iGO Provisional/Laboratory Screw - 1pc
31484	iGO Provisional/Laboratory Wide Screw, Multi-Unit - 1pc
31485	iGO Provisional/Laboratory Wide Screw, Multi-Unit - 5pcs
31486	Multi-Unit/Access iGO Prosthetic Wide Screw - 1pc
31487	Multi-Unit/Access iGO Prosthetic Wide Screw, - 5pcs
31229	Multi-Unit/Access Neo Prosthetic Screw - 1pc
31257	Multi-Unit/Access Neo Provisional/Laboratory Screw - 5pcs
31265	Multi-Unit/Access Neo Laboratory Screw Long - 5pcs
31480	Neo Provisional/Laboratory Wide Screw, Multi-Unit - 1pc
31481	Neo Provisional/Laboratory Wide Screw, Multi-Unit - 5pcs
31482	Multi-Unit/Access Neo Prosthetic Wide Screw - 1pc
31483	Multi-Unit/Access Neo Prosthetic Wide Screw - 5pcs

The *measurement of implant stability* begins with NeoTell™



Resonance Frequency Analysis (RFA)

NeoTell™

Art. No. Description

91060	NeoTell™ Instrument <i>Includes NeoTell™ Instrument, SmartPeg Mount, TestPeg and USB cable, type A – C</i>
91070	SmartPeg Type 18 – Neoss SP (5pcs)
91071	SmartPeg Type 38 – Neoss NP (5pcs)
91072	SmartPeg Type 52 – Neoss Access Abutment (5pcs)
91073	SmartPeg Type 98 – Neoss Multi Unit Abutment (5pcs)
91075	Neoss SmartPeg Mount (5pcs)

New opportunities for your *EasyDigital* workflow



NeoScan

Intraoral Scanners

NeoScan™ 1000

Art. No. Description

91000	NeoScan™ 1000 Intraoral Scanner <i>Includes Handpiece, 4 tips (2×S, 2×L), Holder Set, Manuals, Suitcase. 1 Year warranty included. Extended warranties are available for purchase</i>
91012	NeoScan™ 1000 – Tip L TP102N (2 pcs)
91013	NeoScan™ 1000 – Tip S TP103N (2 pcs)
91024	NeoScan™ 1000 – Warranty extension – 1 year
91026	NeoScan™ 1000 – Warranty extension – 2 years
91027	NeoScan™ 1000 – Warranty extension – 3 years



NeoScan™ 2000

NeoScan™ 2000

Art. No. Description

91001	NeoScan™ 2000 Intraoral Scanner <i>Includes Handpiece, 4 tips (2×S, 2×L), Charging holder, 2 batteries, Stand, Manuals, Wi-Fi adaptor, Suitcase. 1 Year warranty included. Extended warranties are available for purchase</i>
91014	NeoScan™ 2000 – Tip L TP102N (2 pcs)
91015	NeoScan™ 2000 – Tip S TP103N (2 pcs)
91040	NeoScan™ 2000 – Battery (2 pcs)
91032	NeoScan™ 2000 – Warranty extension – 1 year
91033	NeoScan™ 2000 – Warranty extension – 2 years
91034	NeoScan™ 2000 – Warranty extension – 3 years

Digital Impressions

Esthetic Healing Abutments with ScanPeg

All Esthetic Healing Abutments include Titanium Screw and a ScanPeg

Art. No. Description

31360	Esthetic Healing Abutment with ScanPeg – Wide Incisor
31361	Esthetic Healing Abutment with ScanPeg – Narrow Incisor
31362	Esthetic Healing Abutment with ScanPeg – Canine
31363	Esthetic Healing Abutment with ScanPeg – Premolar
31364	Esthetic Healing Abutment with ScanPeg – Molar



ScanPeg

Esthetic Healing Abutments with ScanPeg can also be utilised for digital impressions

Art. No. Description

31365	ScanPeg – 5 pcs Used only with Esthetic Healing Abutments
-------	--



Scan Bodies **NEW**

Art. No. Description

31451	Scan Body Ti Neoss® NP – 1 pc
31452	Scan Body Ti Neoss® SP – 1 pc
31453	Scan Body Ti Neoss® Access – 1 pc
31458	Scan Body Ti Neoss® Multi-Unit – 1 pc
31459	Scan Body Ti Neoss® NP – 5 pcs
31460	Scan Body Ti Neoss® SP – 5 pcs
31461	Scan Body Ti Neoss® Access – 5 pcs
31462	Scan Body Ti Neoss® Multi-Unit – 5 pcs



ScanKey **NEW**


The ScanKey Kits are available as STL files for lab printing. Download the ScanKey at neoss.com.

Art. No. Description

ScanKeyKit
The ScanKey Kit works is compatible with the Scan Body 31452, 31453 and 31458.

ScanKey, NP

The ScanKey Kit – NP is compatible with the Scan Body 31451.



Digital Analogs **NEW**


Art. No. Description

31470	Digital Analog Neoss® NP – 1 pc
31471	Digital Analog Neoss® SP – 1 pc
31472	Digital Analog Neoss® Access – 1 pc
31473	Digital Analog Neoss® Multi-Unit – 1 pc
31474	Digital Analog Neoss® NP – 5 pcs
31475	Digital Analog Neoss® SP – 5 pcs
31476	Digital Analog Neoss® Access – 5 pcs
31477	Digital Analog Neoss® Multi-Unit – 5 pcs
31478	Digital Analog Neoss® Fixation Screw – 1 pc
31479	Digital Analog Neoss® Fixation Screw – 5 pcs

Digital Analog Neoss® Fixation Tool

The Fixation Tool is available as an STL file for lab printing. Download the Tool STL file at neoss.com.





Digital Prosthetics

NeoBase® Abutments SSC (Straight Screw Channel)

also referred to as TiBase

All NeoBase® Abutments SSC include Neo Abutment Screw SP and Neo Laboratory Screw SP



Art. No.	Description	A	B	C
31393	NeoBase® Mono SSC G0.3mm – H3.6mm	3.9mm	0.3mm	4.5mm
31394	NeoBase® Mono SSC G1.5mm – H3.6mm	5.1mm	1.5mm	
31395	NeoBase® Mono SSC G0.3mm – H5.6mm	5.9mm	0.3mm	
31396	NeoBase® Mono SSC G1.5mm – H5.6mm	7.1mm	1.5mm	
31397	NeoBase® Multi SSC G0.3mm – H3.6mm	3.9mm	0.3mm	3.8mm
31398	NeoBase® Mono SSC G0.3mm – H3.6mm, NP	3.9mm		
31399	NeoBase® Mono SSC G0.3mm – H5.6mm, NP	5.9mm	0.3mm	
31400	NeoBase® Multi SSC G0.3mm – H3.6mm, NP	3.9mm		

Include Neo Abutment Screw NP and Neo Laboratory Screw NP

NeoBase® Abutments ASC (Angulated Screw Channel)

also referred to as TiBase

All NeoBase® Abutments ASC include iGO Abutment Screw and iGO Laboratory Screw



Art. No.	Description	A	B	C
31367	NeoBase® Mono ASC G0.3mm – H3.6mm	3.9mm	0.3mm	4.5mm
31368	NeoBase® Mono ASC G1.5mm – H3.6mm	5.1mm	1.5mm	
31369	NeoBase® Mono ASC G0.3mm – H5.6mm	5.9mm	0.3mm	
31370	NeoBase® Mono ASC G1.5mm – H5.6mm	7.1mm	1.5mm	
31371	NeoBase® Multi ASC G0.3mm – H3.6mm	3.9mm	0.3mm	3.8mm
31372	NeoBase® Mono ASC G0.3mm – H3.6mm, NP	3.9mm	0.3mm	
31373	NeoBase® Mono ASC G0.3mm – H5.6mm, NP	5.9mm		
31374	NeoBase® Multi ASC G0.3mm – H3.6mm, NP	3.9mm		


Include iGO Abutment Screw NP and iGO Laboratory Screw NP

Neoss TiBase Abutments (compatible with Dentsply Sirona®)

All TiBase Abutments include Neo Abutment Screw and Laboratory/Titanium Screw



Art. No.	Description	A	B	C
31329	Neoss TiBase N (NB B 3.4 L)	5.3 mm	0.6 mm	4.5 mm
31330	Neoss TiBase W (NB B 4.1 L)		0.6 mm	5.0 mm
31345	Neoss TiBase, NP (FX 3.4 S)	5 mm	0.3 mm	3.9 mm
31331	Neoss ScanPost L (TiBase) Includes ScanPost with Screw SP			
31346	Neoss ScanPost S (TiBase, NP) Includes ScanPost with Screw NP			



Ti NeoLink®

Titanium NeoLinks include Neo Abutment Screw, Neo Laboratory Screw and Straight Copings

Art. No.	Description	A	B	C
31133	Ti NeoLink® Mono	3.9mm	0.3mm	4.5mm
31134	Ti NeoLink® Multi			
31234	Ti NeoLink® Mono, NP			3.8mm

Includes Neo Abutment Screw Ti NP, Neo Laboratory Screw NP and a set of Preformed Plastic Copings



Digital libraries can be found under <https://www.neoss.com/en-us/cad-libraries>

3shape  exocad

Prosthetic Solution Standard and Narrow Platform, overview

Overview Prosthetic Solutions, SP and NP (Ø3.25 – 6.5 Implants)



		Units		Retention		Implant platform		Level			Material					Angulated screw channel	Esthetiline
		Mono	Multi	Cement retained	Screw retained	SP	NP	Implant level	Access level	Multi-Unit Abutment level	Ti	ZrO	Gold	CoCr	Polymer		
Provisional Restorations	Provisional Ti Abutment	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑						
	Esthetic Tissue Former	☑		☑	☑	☑		☑							☑		☑
	Healing Cap									☑							
Final Restorations	Prefabricated	Prepable Ti Abutment	☑		☑	☑	☑	☑			☑						☑
		Zirconia Abutment	☑		☑	☑		☑				☑					☑
		Access Abutment (0 – 30°)		☑		☑		☑			☑						
		Multi-Unit Abutment		☑		☑		☑			☑						
	Digital	NeoBase® (SSC, ASC*)	☑	☑	☑	☑	☑	☑		☑	☑					☑	☑
		TiBase**	☑		☑	☑	☑	☑			☑						
		Ti NeoLink®	☑	☑	☑	☑	☑	☑	☑		☑						☑
		Blanks	☑		☑	☑	☑	☑			☑						☑
		Individual Prosthetics* central milling	☑	☑	☑	☑	☑	☑	☑	☑	☑	☑		☑		☑	☑
	Casting	Gold NeoLink®	☑	☑		☑	☑	☑	☑				☑				☑
		Burnout Abutment	☑	☑	☑	☑	☑	☑	☑								
		CoCr Abutment	☑	☑	☑	☑		☑						☑			
Overdenture	Equator Abutment		☑			☑	☑	☑	☑								
	LOCATOR® Abutment		☑			☑		☑	☑								

* Utilize iGO Abutment Screw and Screwdriver for angulated screw channels, optional for straight screw channels on Individual Prosthetics milled by Neoss

** Compatible with Dentsply Sirona®

Torque and Speed Recommendations, Overview

Neoss Implant System Torque Recommendations (Ncm)

Implants	Screw Taps	Healing			Provisional		Impression	Final restoration					
		Cover Screw	Healing Abutment PEEK and Esthetic Healing Abutment	Healing Abutment Ti	Esthetic Tissue Former	Provisional Ti Abutment	Impression Coping Screw	Neo Abutment Screw	Multi-Unit and Access Abutment Straight & Angulated	Multi-Unit and Neo Access Prosthetic Screw	Wide Screw	Equator and Locator Abutment	
												Implant level	Access level
45 Max	40 – 45	10 Max	10 Max	10 Max	20	20	10 Max	32*	32	20	20	32	20

* 35 Ncm optional for high-load cases

Final restoration with iGO Screws		
iGO Abutment Screw, SP	iGO Abutment Screw, NP	iGO Prosthetic Screw, Access and MUA
32	20	20

Neoss Implant System Drilling/Insertion Speed Recommendations (rpm)

Drills	Screw Taps	Countersinks	Implants	Bone Mill
800 – 2000**	20	800	20 Max	40 Max

** lower speed for larger drills

Symbols, names and abbreviations



SP = Standard Platform Ø3.5 – 6.5 Implants

NP = Narrow Platform Ø3.25 Implants

iGO = utilizes iGO Screw and Screwdriver

Neo = Neoss standard Screw/Screwdriver

SSC = Straight Screw Channel

ASC = Angulated Screw Channel

You can now buy your products online!

Products can also be ordered online at **neoss.com**.

Some products may not be regulatory or commercially available for all markets and available in limited quantities.
Contact your local representative for availability and prices.

Disclaimer of Liability

Neoss products may only be used according to the manufacturers' instructions and recommendations.

The user of Neoss products should determine their suitability for particular patients and indications.

Neoss Limited disclaims any liability, expressed or implied, and shall have no responsibility for any direct, indirect, punitive or other damages arising out of or in connection with any errors in professional judgment or practice in the use or placement of the Neoss products.

Product images are not to scale and are subject to change.

All dimensions are for identification purposes only.

Copyrights, design rights and trademarks:

Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorization of Neoss Limited.

Neoss, 'N' logo, NeoGenix, NeoGen, NeoLink, NeoLoc, NeoBase and Neoss ProActive are registered trademarks of Neoss Group. NeoScan is a trademark of Neoss Group.



Do you want to *lower your inventory costs* and minimize confusion?

NeossONE® is a solution unique to the Neoss implant system – one prosthetic platform, across three implant designs. Simply put, the same prosthetic components fit every implant. The benefits are clear: reducing patient treatment time, optimizing inventory control and minimized product confusion.

- ✓ Reduce steps, saves time
- ✓ Less inventory and maximum flexibility
- ✓ Minimize confusion and mistakes
- ✓ Cost-effective solution for all



NeossONE®



Order your products
24/7 in our Online Store
at shop-us.neoss.com



UNITED STATES

Neoss Inc.

1900 West Park Drive, Suite 165 Westborough, MA 01581

T: +1 866 626 -3677

E: marketingusa@neoss.com

E: NeoScan: neoscan.support.us@neoss.com

For our Terms & Conditions consult your local website on www.neoss.com

Contact your local representative for product availability details.
All products are not available in every market.

10502_25 US 2025-06 © Neoss Limited, 2025. Copyrights, design rights and trademarks, Neoss documents, software and designs may not be reprinted, copied or published in whole or part, without the written authorisation of Neoss Limited.